


Highly sensitive transcriptomic-based pooled CRISPR screening enabled by Spatial Molecular Imager (SMI)

Matthew L. Freedman^{1,2}, Sandor Spisak^{1,2}, Ji-Heui Seo^{1,2}, Talal El Zarif^{1,2}, Shanshan He³, Joe Phan³, Isabel Lee³, Justin Jenkins³, Erin Piazza³, Dae Kim³, Joseph Beechem³
 1-Department of Medical Oncology and 2-Center for Functional Cancer Epigenetics, Dana-Farber Cancer Institute, Boston MA 02215; 3-Nanostring Technologies, Seattle WA 98109



Summary

Spatial molecular imager (SMI) uses nucleic acid hybridization cycles of fluorescent molecular barcodes to enable in-situ measurement of biological targets on intact sample with subcellular resolution. Using SMI, we demonstrate an imaging based pooled CRISPR screens by simultaneously visualizing guide RNA (gRNA) and accurately measuring transcripts of interest within the same cell.

Main Objectives and Method:

- Designed flexible and robust SMI assay capable of simultaneously detecting gRNAs and quantifying RNA expression in the same cell that can scale in target plex (up to 1000) and number of cells interrogated (>100K)
- Using LNCaP prostate cancer cell line, applied epigenetic perturbations to an androgen receptor (AR) enhancer to validate impact on the AR gene and its downstream targets (55 gRNAs + 37 genes + controls)

Key Results:

- High Resolution** – Unambiguous visualization and digital detection of gRNAs in single cells
- High Sensitivity** – Demonstrated by:
 - Detection of key AR related target genes that are typically low expressor (TPM ~ 1-10)
 - Quantitation of AR related target genes on over 90% of the cells analyzed (~10% drop-out rate)
 - High detection sensitivity of rare modified cells (below 1%) was observed
- Good concordance** – Gene expression changes measured by SMI vs. Z-scores from RT-PCR drop-out assays


Key advantages:

- High plex – Flexible panel design to visualize and measure up to 1000 gRNAs and genes within single cells
- Accurate quantitation – Decreased 3'-end transcript bias and use of both endogenous and exogenous controls
- High throughput with simple workflow – Tens to hundreds of thousands of cells analyzed per SMI run


For research use only. Not for use in diagnostic procedures.

Introduction: Spatial Molecular Imager

SMI is a single instrument solution for subcellular spatial analysis



SMI's cyclic readout enables high-plex detection




x8 Enabling highest sensitivity


Barcode encoding scheme: 64 bits (4-colors x 16 cycles)

* More technical details on SMI on Beechem et al. AGBT 2021 poster


Fundamental Challenge – Connecting regulatory elements to their genes




Method Overview: Pooled CRISPR screening using transcriptomic readout on SMI




SMI produces high data density in virtually every cell




Proof-of-concept: Measuring effect of genomic perturbation by RNA readout using SMI



SMI data: Statistical analysis of transcript distribution between two gRNA designated cells [cells detected with HPRT1 gRNA vs. non-target gRNA (CTRL)]



Visualize and quantify genomic perturbation at single cell resolution on tens of thousands of cells using SMI



Conclusions

Using SMI technology, we demonstrated the ability to conduct highly multiplexed gRNA pooled screening, which seamlessly allows simultaneous identification of gRNA and accurate digital measurement of transcripts of interest at single cell resolution. Concordance between this approach with existing RT-PCR based drop out assay validates the ability to perform CRISPR screens using gene expression as a readout. This approach, with further refinement would be transformative for understanding regulatory element biology because it allows a more direct and accurate way to functionally connect regulatory element-to-target gene relationships.

