Spatially resolving the host immune-epithelial landscape in lung tissue from patients infected with SARS-CoV-2 with proteomic and whole transcriptomic profiling nanoString

(4) Dispense cliscs into 56 nuel plote

Bronchiolar Tissu

Gene of Interest

UMAP 1

Disease:Healthy LFC of (A) anti-viral gene responses and (B)

markers of inflammation in ARDS (blue) and COVID-19 (red)

plot illustrating the unique global response between healthy.

COVID-19, and ARDS alveolar lung tissue

alveolar and bronchiolar tissue. (C) UMAP dimension reduction

Figure 3:

CXCL10 CXCL8 STAT1 STAT2 SFTPD MUC1

Region

Patient

Aveciar

Bronchiola

Disorganiz

ARDS-2 440

40021441

COMD_050

COMD 073

AFFK138

6 Construct Liketry

Jamie L Everman¹, Michael T Montgomery², Alison VanSchoiack³, Tyler Hether³, Michael Leon³, Kristina Sorg³, Wes Goldman³, Yan Liang³, Sarah E Warren³, William J Janssen⁴, Max A Seibold^{1,5} 1. Center for Genes, Environment, and Health, National Jewish Health, Denver, CO, USA. 2. Department of NanoString Technologies Genetics, Stanford University, Palo Alto, CA, USA. 3. NanoString Technologies, Inc., Seattle, Washington, USA. 4. Department of Medicine, National Jewish Health, Denver, CO, USA. 5. Department of Pediatrics, National Jewish Health, Denver, CO, USA. 530 Fairview Avenue North, Seattle, WA 98109



(PanCK) and immune cells (CD45), or for nCoV2019-S by

(B) Sankey plot illustration of the distribution of selected

ROIs for DSP analysis across health status and anatomical

RNA-FISH (RNAscope™), nuclei are labeled with DAPI.

region of the lung tissue

Inflammation- and disease-specific immune cell infiltration in COVID-19

- · A distinct inflammatory, immune cell, and epithelial cell response can be seen between mild/moderate and severe inflammatory states of COVID-19 (gene expression; panel A)
- · COVID-19 drives a unique tissue-associated cell type profile in ROIs with mild/moderate or severe inflammation compared to ROIs from healthy subjects (enriched cell types; panel B)



· Proteomic profiling demonstrates T-cell, B-cell, and fibroblast protein signals that are higher in expression in COVID-19 alveolar tissue compared to healthy or ARDS alveolar lung tissue



Figure 4: (A) Gene expression changes or (B) lung and immune cell type enrichments plotted across healthy and COVID-19 alveolar lung tissue based on mild/moderate versus severe inflammation: color = p<0.01: (C) Heatmap of expression levels of select assayed protein targets in alveolar tissue from COVID-19 (n=2), healthy (n=1), and ARDS (n=1) subjects.

Conclusions

- · COVID-19 drives a robust anti-viral gene expression response in both alveolar and bronchiolar lung tissue Lung tissue from COVID-19 patients demonstrates a unique inflammatory expression profile, in both cytokine
- signals and inflammatory pathways, compared to lung tissue from ARDS patients
- Differing states of COVID-19-driven inflammation alters immune cell type presence in COVID-19 tissue
- Inter-donor and intra-donor heterogeneity in COVID-19 subjects is significant and drives differing levels of inflammation and anti-viral responses, immune cell infiltration, and gene and protein expression patterns in both alveolar and bronchiolar airways and lung tissue

Acknowledgements/Disclaimers

We would like to thank National Jewish Health, Mt. Sinai Hospital, University of Oklahoma, and University of California San Diego for providing healthy, ARDS, and National Jewisl COVID-19 lung tissue sections for this study. Health NanoString products are For Research Use Only and not for Science Transforming Life use in diagnostic procedures.



FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. www.nanostring.com | info@nanostring.com © 2021 NanoString Technologies, Inc.

3 (12)

CVD19