


## Summary

In this study, we used the GeoMx<sup>®</sup> Whole Transcriptome Atlas and RNAscope to visualize repeat elements and characterize the impact of repeat elements-presenting regions of pancreatic cancer. Then, we used CosMx<sup>®</sup> SMI to collect gene expressions of spatially-annotated single cells in both repeat elements-positive and repeat elements-negative regions. SMI uses fluorescent molecular barcodes to enable in-situ measurement of biological targets on an intact tissue sample, such as formalin-fixed paraffin-embedded tissues using a 1000+plex RNA panel. Using the RNA expression data from SMI, we were able to spatially visualize how EVs encapsulating repeat elements impact the surrounding CAFs and other cells in the microenvironment. Using cell type deconvolution algorithms, we characterize the changes in gene expression within cell populations using the GeoMx as well as individual cells using SMI. Using nearest neighbor analyses, we are able to characterize expression profile changes between CAFs and immune cells that are close and distant from repeat expression tumor cells.

These results have provided a new lens in viewing the impact of repeat element tumor expression and mechanistic role in the microenvironment via EVs, which have the opportunity to identify novel therapeutic avenues and biomarkers in pancreatic cancer.


## Introduction



Repeat sequences have been shown to play a role in wide range of diseases including cancer, chronic viral infection, and neurodegeneration. A recent Phase II clinical trial of pancreatic cancer indicating that expression of certain repeat RNAs are correlated with response to combined immune checkpoint blockade and sublethal radiation therapy. More mechanistic studies have now revealed that these repeat elements are transcribed and often packaged in extracellular vesicles (EVs) which can be released from the tumor cells and interact with cells in the tumor microenvironment including cancer-associated fibroblasts (CAFs). Thus, it is important to understand the spatial distribution of repeat RNAs found in tumor cells and the surrounding CAFs that can provide mechanistic insight into their consequence in shaping the cellular response to tumors. Of particular interest is the impact they have on CAF subtypes that have been defined by single cell RNA-seq (scRNA-seq), which include the inflammatory (iCAF) or myofibroblastic (myCAF) phenotypes. Given the limitations of bulk RNA-seq and scRNA-seq in understanding the precise spatial distribution of this phenomenon, we used the NanoString GeoMx and Spatial Molecular Imager (SMI) to visualize these repeat elements and study their impact on the surrounding cells.


## Experimental design

Generate multi-omics data sets to investigate roles of repeat RNAs in PDACs




## CosMx SMI assay overview

CosMx SMI is a single instrument solution for subcellular spatial analysis




## High pleX CosMx SMI panel



## CosMx SMI visualizes 965 gene targets in Pancreatic Cancer tissues


CosMx SMI spatially resolve 965 target genes across 3 PDAC tissues and 1 normal pancreatic tissue.




4 PDAC tissues	
Total number of cell detected	337,969
% cells (>20 RNA counts, %)	98
Total number of transcripts detected (count)	114,966,012
% transcript assigned to cell (%)	98
Total FOV number	98
Mean transcript/cell (count)	263
Mean transcript/cell/target (count)	263
Mean neg probe/cell/target (count)	0.85

## Customized CosMx cell typing analysis identifies 21 cell types of PDAC tissues and visualize their spatial locations


Gene expression signatures specific to cell types are defined using scRNA-seq data. Nanostring cell typing algorithm determines cell types matching to a tissue type. CosMx SMI retains x,y coordinates of cells unlike dissociated scRNA-seq data.





## Repeat RNA-expressing cancer cells are losing keratin expressions, suggesting that repeat RNA promotes EMT





D10 (PDAC) RNA expression and cell types  
 RNA: LINE1\_ORF2, KRT19  
 Cell: Cancer cell Undifferentiated Epithelial-like cancer cell



## GeoMx DSP profiles whole transcriptome expressions in different tissue compartments




## Contour-1-2 has the highest repeat RNA expressions



## High-repeat-RNA tumor-microenvironment regions have upregulated genes which are associated with cytokine and inflammatory signaling

- Pathway Names: Cytokine Signaling in Immune system, ECM proteoglycans, Extracellular matrix organization, Interferon gamma signaling, Interleukin-4 and Interleukin-13 signaling, KRT, Non-integrin membrane-ECM interactions
- Pathway Names: Cytokine Signaling in Immune system, Interferon gamma signaling, Interleukin-4 and Interleukin-13 signaling



## Conclusions

- CosMx SMI is a single instrument solution for subcellular spatial analysis and provides sub-cellular resolution of 1000+plex transcriptomic information.
- Using the RNA expression data from CosMx SMI and GeoMx DSP, we were able to spatially visualize how repeat elements impact cancer cells and their potential delivery via EVs in affecting surrounding cells in the microenvironment.
- Using pathway analysis, we characterize the changes in gene expression within cell populations using the GeoMx and provide single cell resolution validation using CosMx SMI.
- Cellular subsets with gene expression profiles defined using single-cell sequencing-based multiomics can be spatially resolved using SMI to uncover novel cell-cell interactions.
- These results have provided a new lens in viewing the impact of repeat element tumor expression and mechanistic role in the microenvironment via EVs, which have the opportunity to identify novel therapeutic avenues and biomarkers in pancreatic cancer.

## Reference

Elyada et al., Cancer Discovery (2019)  
 Sodergren et al., J Cancer Res Clin Oncol (2020)

