

# #202 Integrating single cell and spatial gene expression profiling of mouse organogenesis to identify and localize unknown cell types

Chengxiang Qiu<sup>1</sup>, Eva Nichols<sup>1</sup>, Stephanie Zimmerman<sup>2</sup>, Alexa Lasley<sup>2</sup>, Margaret Hoang<sup>2</sup>, Joe Beechem<sup>2</sup>, Jay Shendure<sup>1</sup>  
1. Department of Genome Sciences, University of Washington, Seattle, WA, USA, 2. NanoString® Technologies, Seattle, WA, USA








**Abstract**

Mammalian organogenesis is a remarkable process, whereby cells within the post-gastrulation embryo continue to rapidly proliferate while giving rise to the diverse cell types of each organ system, specified by molecular programs that are precisely regulated in time and space. Single cell RNA-sequencing of whole embryos during mouse embryogenesis and organogenesis is yielding unprecedented detailed views of mammalian development, for example revealing hundreds of unique cell types defined by gene expression. Although many methods exist for identification of cell types defined by scRNA-seq, annotating cells remains a manual and challenging process. In this work, we sought to leverage spatial gene expression data of mouse organogenesis to validate annotations and localize uncertain cell populations to specific tissues or regions.




We used high-resolution scRNA-seq data from mouse development -- specifically, single nucleus transcriptional profiling of millions of cells done by 3-level combinatorial indexing. This "whole organism profiling" was performed on staged embryos in 2 to 6 hour increments from gastrulation to birth. The resulting single cell profiles were processed by conventional methods and, although an initial round of manual annotation based on marker genes and earlier generations of atlases was fruitful, many ambiguities remained. To address these in part, we integrated matched timepoints with spatial whole transcriptome profiles of specific anatomical structures of four timepoints of mid-gestation mouse development generated using the GeoMx® Digital Spatial Profiler (DSP). We used a cell type deconvolution algorithm to estimate the abundance of each cell type in each region profiled by DSP and validated that known cell types such as tissue-specific epithelial cell subtypes localize to the correct anatomical structures with high accuracy. We then used this method to further map the cell trajectories derived from the lateral plate mesoderm, populations which have limited research and are therefore challenging to annotate.

Next, we applied this method to understand how dysregulated cell lineage contributes to organ malformation in a developmental mutant. Absence of embryonic macrophages due to colony stimulating factor 1 receptor (CSF1R) deficiency causes bone and brain deformities as well as perinatal lethality in mouse and humans, suggesting important functions in organ formation. We performed massively scalable RNA single-cell transcriptomics (via single-cell combinatorial indexing RNA sequencing) and GeoMx DSP on E18.5 embryonic tissue sections of wildtype and CSF1R-deficient mutant littermates. We find differential cell type abundance in both the scRNAseq and between matched spatial regions in wild type and mutant in a wide variety of tissues, suggesting that organs beyond bone and brain are impacted by embryonic macrophage loss. In conclusion, this work provides a framework for integrating spatial data with scRNAseq in an automated pipeline to add spatial annotations to unknown cell types in normal and pathological samples.


## Methods: single cell RNAseq and spatial whole transcriptome profiling of mouse organogenesis



## Application to a developmental mutant: sciRNAseq and spatial profiling of *Csf1* deficient embryos



## Mapping major cell types spatially and identifying tissue-specific differential cell abundance in *Csf1* mutants



**Conclusions**

- High resolution sciRNAseq identifies novel cell types during dynamic processes such as embryonic development
- Spatial whole transcriptome profiling with GeoMx WTA can localize and annotate novel cell types identified by scRNAseq
- In developmental mutants, single cell RNAseq plus spatial gene expression data can link observed pathology to changes in cell abundance, cell localization, or cell state

**References**

- Cao, J. et al. *Nature* 566: 496–502 (2019).
- Martin, B. K. et al. *Nat Protoc* 18: 188–207 (2023).
- Zimmerman, S.M. et al. *Genome Res* 32(10): 1892–1905 (2022).
- Denaher P. et al., *Nat Comm* 13(1): 385 (2022)
- Chen, A. et al. *Cell* 185(10): 1777–1792 (2022)
- Hume, D. A. et al. *J Leukoc Biol* 107: 205–219 (2020).

Scan here to download or learn more