127 Large scale and high-throughput spatial-omics to deeply characterize cancer and normal tissue

Sarah E. Church, Margaret Hoang, Gokhan Demirkan, Kyla Teplitz, Joachim H. Schmid, Joseph M. Beechem





To understand heterogeneity of human tissue and identify biomarkers for predicting response to therapies it is necessary to evaluate thousands of patient samples. Here we present a process for high-throughput spatial analysis of the whole transcriptome and over 100 proteins for thousands of clinically relevant cancer and normal samples using the GeoMx[®] digital spatial profiling (DSP) platform. From these analyses we are able to characterize the tumor microenvironment in separate compartments of the interepithelial, stromal and immune cell rich areas of multiple solid tumor indications and subtypes. Similarly, in normal tissue we are able to deeply characterize tissue specific histology such as glomeruli in kidney, hippocampus in brain or peyer's patches in colon tissue.

Abstract

Here we present a method to spatially analyze 96 FFPE samples on a DSP instrument per week using established automated procedures for each step of the DSP workflow. This process incorporates region of interest (ROI) and area of illumination (AOI) selection, slide processing, the DSP instrument run, sequencing preparation, sequencing and data analysis. Pre-selection of up to 12 ROIs/AOIs per sample is performed on sequential H&E slides, which then can be rapidly overlayed onto the DSP immunofluorescent scans, easily placed and automatically segmented on specific tissue such at cytokeratin-positive tumor and cytokeratin-negative stroma. Slide processing is done automatically using the Leica BOND system for 32 slides with 3 samples per slide. All annotations chain of custody can be tracked during the process using a barcode reader. After DSP on-instrument processing, collection plates can then be stored until ready for automated pre-sequencing processing. Following sequencing, data is processed and uploaded to AtoMx[™] Spatial Informatics Platform for data QC, analysis, collaboration and storage.

This high throughput workflow spatially analyzing both protein and whole transcriptome gene expression allows for the creations of databases and consortiums containing meta-data for large numbers of tissues. These databases can be used to train new AI methods, to deeply evaluate tissue neighborhoods, and in the case of clinical trials, to discover biomarkers to understand mechanism of disease and potentially help improve selection of therapies in the clinic.

32 slides in two batch of 16 slides stained with protein or WTA	NGS RNA readout fully automated from start to finish including fluorescent morphology markers. Semi-automated protein prep	Slides stored and stable for profiling for 2 weeks	



large cohorts of cancer and normal samples that is based on tissue architecture and can be easily combined across multiple instruments, institutions and users to generate spatial multi-omic databases to better understand mechanism of disease

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. www.nanostring.com | info@nanostring.com ©2023 NanoString Technologies, Inc.