^{#615} Spatially resolved expression of T cell receptors elucidates spatial relationships between T cells, immune infiltration, and cancer-associated pathways

Katrina van Raay¹, Jason Reeves¹, Erin Piazza¹, Hargita Kaplan², John Vivian², Francis Fernandez², Margaret Hoang¹, Joseph Beechem¹ ¹NanoString® Technologies, Seattle, WA ²Atreca, Inc., San Carlos, CA *Correspondence: jbeechem@nanostring.com

Abstract

Spatial distribution of T cells is key in understanding the escape of tumors from immune surveillance via the adaptive immune response, including interactions between immune cells and the surrounding tumor microenvironment. T cells are critical to the adaptive immune response to pathogens and cancers, mediating an antigen-specific response through both specificity and diversity of T cell receptor (TCR) clonotypes. Many methods exist to determine specific clonotypes and overall TCR diversity present from bulk tissues or sorted cell populations; however, nearly all fail to capture spatial orientation and arrangement of T cells engaging with their microenvironment, and most require large amounts of starting material from precious samples. Here, we present a TCR expression profiling panel for the GeoMx® Digital Spatial Profiler that can be combined with the GeoMx Cancer Transcriptome Atlas (CTA) or Human Whole Transcriptome Atlas (WTA) on archival formalin-fixed paraffin embedded (FFPE) tissue specimens. This represents the first commercial spatial expression profiling assay for the simultaneous quantification of TCR constant, variable, and joining segments in situ.

We show reliable sensitivity and specificity (>90%) with respect to orthogonal sequencing and robust detection of TCR chains with evidence of clonal expansion and CD8 infiltration across tumor regions in colorectal cancer tissue. These events also corresponded to increased signatures of exhaustion from the T cells and suggest that the T cells resident in or near the tumor are tumor-specific and poised for activation via checkpoint blockade. Signaling pathways and tumor-specific signatures were also evaluated to look for mechanisms through which tumor cells respond to T cell infiltration. We further validated the performance of the TCR probe pool in cell pellet arrays with orthogonal TCR sequencing, tonsil and colorectal cancer tissues.

Together, the combination of our TCR add-on panel with the CTA or WTA illuminates T cell phenotypes, signaling pathways, population dynamics, and transcriptomic changes, yielding an unparalleled view of the T cell response in any context.



alysis is performed with the GeoM 6 Data analysis is performed with the Geo Data Analysis Suite and custom software ollected probes are counted on nCounte

Tagged Oligonucleotide Chemistry

5 Analysis System or by a Next-Gen Sequence

GeoMx Digital Spatial Profiler (DSP) uses oligonucleotides which hybridize to target mRNAs to quantitatively read out DNA tags which are selectively released *in situ* by specifically shining UV light into certain regions of the tissue.

compartments, complex morphological structures, and cell populations dispersed throughout a tissue.



GeoMx DSP enables a *tailored* profiling approach to tackle specific questions



A variety of pre-set tools already exist to create Regions of Interest (ROIs) to your specifications. ROIs may also be hand-drawn, creating limitless possibilities for sampling designs

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. www.nanostring.com | info@nanostring.com

© 2023 NanoString Technologies, Inc.

alpha, beta, gamma, delta chains, and:













nanoString

NanoString Technologi 530 Fairview Avenue North, Seattle, WA



- T-cell receptor gene segments were confirmed to be sensitively and specifically captured by GeoMx DSP
- Diverse T-cell populations detected from melanoma, along with sample-restricted detection of γ/δ T-cells
- GeoMx simultaneously differentiated α/β from γ/δ -derived T cell lymphomas on a tissue microarray and identified genes enriched

