GeoMx[®] Digital Spatial Profiler

Grant Support Package



NanoString Technologies, Inc. approves the use of images, graphs/charts, methods, and descriptions found in this grant package with proper citing for grant proposal purposes. **Please cite:** [Image, data, methods, description] courtesy of NanoString Technologies, Inc.

Section 1: Information to provide in grant application Platform Overview

GeoMx[®] Digital Spatial Profiler (DSP) is the solution for spatial multi-omic profiling.

DSP enables morphology-driven, high-plex profiling on a single FFPE or fresh frozen tissue section. DSP technology combines standard immunofluorescence techniques with oligonucleotide barcoding for readout via NGS or on an nCounter[®] Analysis System. The system specifications are as follows:

- 1. Throughput: >10 sections per day
- 2. Multi-analyte: RNA and Protein
- 3. Profiling via geometric, contour, gridded, and morphology marker-driven segmentation.
- 4. Minimum Segment Size Recommendation:
 - Target detection is dependent on the size and number of cells in an area of interest and the expression level
 - For small cells (e.g., T Cells or B cells): 10-20 cells for protein, 50-200 cells for RNA
 - For large cells (e.g., astrocytes, neurons): 1-5 cells for protein and 5-20 for RNA
- 5. Multiplexing capability:
 - 96 targets with nCounter readout
 - 20,000+ targets with NGS readout
- 6. Quantitative resolution: 6 logs of dynamic range

Technology Overview & Workflow

The GeoMx® Digital Spatial Profiler (DSP) is a novel platform developed by NanoString Technologies. DSP enables spatially resolved, high-plex (10s - 10,000s of targets) digital quantitation of proteins and mRNA in tissue. The assay utilizes unique reagents (antibodies or RNA probes) coupled to UV photocleavable oligonucleotide barcodes. After incubation/hybridization of the GeoMx reagents along with visualization reagents to slide-mounted formalin fixed paraffin-embedded (FFPE) or fresh frozen tissue

GeoMx DSP Workflow

sections, the DSP scans the slides and gives a high quality image for region of interest selection. The oligonucleotide tags are released from user-selected regions of the tissue with focused UV light based on the user region and segmentation selections. The digital micromirror device of the DSP instrument tunes the UV light with 1 micron resolution. This allows great flexibility for region selection including irregularly shaped and noncontiguous segments. Released tags are counted on the nCounter Analysis system or sequenced on an Illumina® next-generation sequencer (NGS). Finally, counts are mapped back to the tissue location by the software, resulting in a spatially-resolved digital profile of protein or mRNA abundance.

- A 5-10µm thick formalin-fixed, paraffin-embedded (FFPE) or fresh frozen tissue section is stained with a cocktail of standard morphology markers (i.e., fluorescent antibodies or RNAScope probes) combined with either oligo-tagged antibodies (for protein detection) or oligo-tagged ISH probes (for RNA detection). Oligo tags are attached to the antibody and ISH probes via a UV-photocleavable linker. Staining can be performed manually or automated using the Leica BOND RX Fully Automated Research Stainer.
- Regions of interest (ROI) on the tissue section are selected by the user using the control software. Using fluorescent morphology markers and our autosegmentation algorithm, users can further divide an ROI into individual segments, such as specific tissue compartments or cell populations.(Please refer to section: GeoMx DSP ROI selection).
- 3. After ROI/segment selection, oligos from the selected regions are released upon focused exposure to UV light.
- 4. The photocleaved oligos are aspirated via a microcapillary tube and stored in microplate wells.
- Lastly, the photocleaved oligos from the spatiallyresolved ROIs in the microplate are (1) quantified using standard NanoString nCounter or next-generation sequencing Illumina[®] workflows.



GeoMx DSP Instrumentation and Assays

Spatial biology experiments with GeoMx require the DSP instrument for slide scanning, ROI selection, and oligo collection. Downstream of the DSP, there are two instrumentation platforms for quantitative readout: The NanoString nCounter[®] platform or a standard library preparation workflow for analysis on an Illumina[®] NGS platform.

The GeoMx DSP platform is commercially available from NanoString Technologies. DSP results correlate with those from immunohistochemistry; high concordance is also observed between DSP results generated using both antibodies or RNA probes. Oligo conjugation does not interfere with antibody or RNA probe specificity and conjugation ensures robust quantification of protein and RNA abundance. Additionally, analytical validation data has demonstrated multiplexing compatibility. Lastly, any commercially-available antibody can be custom-labeled with oligo tags and used in a DSP protocol. (referenceinternal data)



GeoMx DSP ROI Selection

The DSP software enables the user tremendous flexibility in the ROI selection strategy to profile tissue. In particular, highlighted below are five unique profiling modalities designed to interrogate tissue samples.

The segmentation and cell-type specific approaches are achieved through the use of masks in the software based on the fluorescent marker staining of the tissue (e.g. PanCK+ vs. PanCK-). UV light then specifically targets the ROI with multiple illumination events, enabling the separate profiling of complex and interdigitated cell populations within a selected ROI (e.g., tumor vs. stroma).



GeoMx DSP Data Analysis

The GeoMx DSP comes with built-in data analysis software that enables users to easily visualize and interpret their data. The intuitive interface displays high-resolution images side-by-side with quantitative profiling data. The connection between the image and profiling data is preserved allowing users to interact with either dataset in real time. For example, users can select a data point within a data visualization and the software will automatically highlight the ROI or segment from which the expression data was generated. The data analysis workflow includes QC and normalization operations, as well as several data visualization methods for various types of scientific and statistical inquiries.



Users can auto-generate figures such as heatmaps, box plots, and bar graphs, to visualize the distribution of protein or gene expression across multiple tissue compartments (e.g., tumor versus stroma); discover statistical changes in expression between experimental groups (e.g., t-tests, linear mixed models) that can also be visualized in volcano and PCA plots; visualize large data matrices using dendrograms and clustering analyses to discover biomarkers, and perform pathway analyses to understand disease mechanisms of action and identify novel targets.

DSP data analysis is also available in the cloud with the AtoMx[™] Spatial Informatics Platform (SIP). AtoMx SIP enables users to perform image analysis & data visualization, utilize scalable cloud compute and storage, and collaborate globally with researchers across the globe.



GeoMx DSP Consumables

The GeoMx DSP is a multi-analyte platform that can spatially resolve both protein and RNA. Protein and RNA can be separately resolved on two serial tissue sections or resolved simultaneously from a single tissue section with the Spatial Proteogenomics workflow. Validated off-the-shelf assays include content for immuno-oncology, immunology, and neuroscience including the Whole Transcriptome Atlas for human or mouse.

Protein assays are available for either nCounter or NGS readout. Assays are designed for flexibility with a ~20 plex core panel that includes additional positive and negative controls. Additional ~10-plex modules are available and can be mixed and matched with the core based on an individual researcher's needs (up to 96-plex total for assays utilizing nCounter readout and 90+ plex for assays utilizing NGS readout.). Users can also add up to 5 custom targets with a DIY protein barcoding kit or up to 10 custom targets of interest using NanoString's Protein Barcoding Service (more details for custom antibodies below). A Core panel must always be run.

The following content is commercially available for nCounter readout

- <u>Human IO Protein –</u> Immune Cell Profiling Core, IO Drug Target Module, Immune Activation Status Module, Immune Cell Subtyping Module, Pan-Tumor Module, Cell Death Module, AKT/PI3K Signaling Module, and MAPK Signaling Module
- 2. <u>Mouse IO Protein</u> Immune Cell Profiling Core, IO Drug Target Module, Immune Activation Status Module, Immune Cell Subtyping Module, Cell Death Module, and Pan-Tumor Module
- 3. <u>Human Neuroscience Protein Neural Cell Profiling</u> <u>Core</u>, Alzheimer's Pathology Module, Parkinson's Pathology Module, Alzheimer's Pathology Extended Module, Autophagy Module, and Glial Cell Subtyping Module
- 4. <u>Mouse Neuroscience Protein:</u> Neural Cell Profiling Core, Alzheimer's Pathology Module, Parkinson's Pathology Module, Alzheimer's Pathology Extended Module, Autophagy Module, and Glial Cell Subtyping Module

The following content is commercially available for NGS readout:

 <u>Human and Mouse:</u> Protein Core, Immune Cell Subtyping Module, IO Drug Target Module, Immune Activation Status Module, Pan-Tumor Module, Myeloid Module, MAPK Signaling Module, PI3K/AKT Signaling Module, Neural Cell Typing Module, Alzheimer's Pathology Module

All GeoMx DSP protein assays are validated with a preand post-conjugation screen, a sensitivity/specificity screen using FFPE cell pellets, an interaction screen, and a TMA screen using normal and tumor TMAs (see Supporting Data).

Custom antibodies:

Custom Antibodies from Abcam: Abcam provides 1-5 custom antibodies to extend GeoMx protein assays. To add custom antibodies to a GeoMx protein assay, users place the order directly through Abcam, using one of Abcam's IHC-validated recombinant rabbit monoclonal antibodies, and the form found at: <u>www.abcam.com/dsp</u>. Abcam custom antibodies are tested in a pre- and post-conjugation screen.

Custom Antibodies from NanoString: NanoString provides 1-10 custom antibodies through the NanoString Protein Barcoding Service (PBS). Users purchase IHC-validated antibodies directly from a vendor and ship them to NanoString for barcoding. NanoString custom antibodies are tested for the presence of the conjugated tag, but are not tested with IHC or further validated.

A DIY Protein Barcoding Kit can be used instead of purchasing custom barcoded antibodies from Abcam or NanoString and up to 5 antibodies can be added to the DSP protein assays for either nCounter or NGS readout.

Guidance on further validating custom antibodies can be provided upon request.

RNA Content:

 <u>Human and Mouse IO RNA Assay for nCounter</u> <u>Readout</u>- Immune Cell Pathways Panel with 84 targets (including controls) and room to add up to 10 custom targets (+2 additional controls) of interest for a total of up to 96 targets. Fixed panel content covers characterization of the global immune response, microenvironment immune activity, tumor reactivity, and the Tumor Inflammation Signature.

2. Human Cancer Transcriptome Atlas (CTA) for NGS:

1,800+ mRNA gene targets that cover important aspects of the immune response, tumor biology and the microenvironment. The CTA curated content is deeply annotated with 55 major pathways, immune cell subsets, and clinically-derived gene sets such as the Tumor Inflammation Signature (TIS) and PAM50 tumor subtype geneset. Users can add up to 200 custom mRNA gene targets of interest.

- 3. <u>Canine Cancer Atlas (CCA) for NGS Readout:</u> Built with input from experts in canine comparative oncology, the Canine Cancer Atlas contains overlapping content with the nCounter[®] Canine IO Panel, allowing you to seamlessly transition your comparative oncology studies from bulk gene expression analysis to spatial transcriptomics. Profile 1,962 RNA targets with spatial resolution and spike in up to 200 custom RNA targets to add in probes for non-coding RNA, synthetic DNA, or exogenous genes. Targets include those involved in the onset, progression, and immune response of common canine cancer such as melanoma, osteosarcoma, lymphoma, urothelial carcinoma, and glioblastoma.
- 4. <u>Human Whole Transcriptome Atlas (HuWTA) for NGS</u> <u>Readout</u>: 21,000+ protein encoding mRNA targets encompass the whole human transcriptome and is relevant for any type of research. The HuWTA takes a targeted approach to quantify transcripts in every biological pathway and across all gene expression levels. It has been pruned of uninformative high expressors like ribosomal proteins from the target list, so that reads are comprised of informative protein coding genes. Users can add up to 200 custom targets of interest to count unique transcript variants such as isoforms, long noncoding RNA and exogenous transcripts.
- 5. Mouse Whole Transcriptome Atlas (MuWTA) for NGS <u>Readout</u>: 18,000+ protein encoding mRNA targets encompasses the whole mouse transcriptome and is relevant for any type of research. The MuWTA takes a targeted approach to quantify transcripts in every biological pathway and across all gene expression levels. The probes have been designed toward BalbC, C57B6, and NOD strains and has considerable cross-reactivity with other strains. It has been pruned of uninformative high expressors like ribosomal proteins from the target list, so that reads are comprised of informative protein coding genes. Users can add up to 200 custom targets of interest to count unique transcript variants such as isoforms, long non-coding RNA and exogenous transcripts.

All NanoString RNA probes are validated in a bioinformatic screen and undergo technical QC (see Supporting Data).

Training and support

Upon purchase and acquisition of the GeoMx DSP instrument, NanoString will provide on-site installation and calibration of the instrument by a trained field engineer. Once the instrument has been qualified, a Product Application Scientist (PAS) will be on site to provide a comprehensive introduction to the GeoMx DSP platform. The introduction will cover routine use of the instrument, experimental design, project consultations, and data analysis. In addition, the PAS will train the lab scientist to ensure that they have hands-on experience in operating the instrument, implementing the selected GeoMx protocol, and running an experiment from end-to-end. Lastly, the PAS will provide data analysis training with the GeoMx Data Analysis software or the AtoMx SIP so that the customer is familiar with all data analysis methods to process, analyze and interpret their data. Upon successful completion of the training, the customer can access support through the NanoString GeoMx Customer Experience Team regarding experimental design, instrumentation, consumables, and software.



Tips to maximize success with the grant application process

Emphasize Collaboration

The GeoMx DSP software and AtoMx SIP enables collaboration and secured remote-access. The GeoMx instrument makes it possible to securely access your data from anywhere and at anytime through a web browser. The software also provides administrators the ability to assign roles and permissions to their team members and collaborators. The platform allows multiple users within large labs, core facilities, or contract research organizations to securely manage and share their data.

Enabling the DSP workflow in your research environment will facilitate deep and meaningful collaborations between scientists, pathologists, and biostatisticians. Hypothesis-driven study design, pathology-driven ROI selection, and data-driven analyses will enable translational discoveries and contributions to the fields of science and medicine. Acknowledging these collaborative opportunities and providing details on how the working relationships will be fostered can increase the success of your grant application.

Benefits to local research environment

The digital spatial profiler (DSP) quantifies biology with spatial context. Any tissue that can be sectioned including Formalin-Fixed Paraffin Embedded (FFPE) and fresh frozen tissues can benefit from the platform. Our reagents are designed towards human and mouse, but our reagents have appreciable and useful cross-reactivity with non-human primate and rat respectively. Information on cross-reactivity is available upon request. The DSP is relevant to any individual using tissue slides to study biology, making the technology a valuable addition to the local research ecosystem.

Publications

The utility of the GeoMx DSP system is demonstrated already in over 200 peer-reviewed publications, many in leading journals. Please click <u>here</u> to view current GeoMx DSP Publications.

Supporting Data

Reproducibility of GeoMx DSP

All GeoMx antibodies and RNA Probes go through a robust validation process to ensure they will perform as expected and researchers can use GeoMx reagents off the shelf without having to spend time validating panels. Additionally, serial sections analyzed on the GeoMx DSP for protein (Figure 1) or RNA (Figure 2) expression demonstrate high reproducibility. For additional reproducibility data see: Gupta et al., 2020 (Rimm); Merritt et al., 2020



FIGURE 1: Two serial tonsil sections demonstrate high concordance with 24 co-registered ROIs



FIGURE 2: GeoMx® DSP High-Plex RNA was performed on 5 µm FFPE sections. High reproducibility was demonstrated between experimental replicates with serial sections with a pre-commercial panel of 4,998 probes.

Ordering Information

Ancillary Products

Product	Description	Catalog Number
GeoMx* Digital Spatial Profiler Slide Tray	Additional 4-slide holder	GMX-DSP-TRAY

GeoMx Digital Spatial Profiler

Product	Description	Catalog Number
GeoMx Digital Spatial Profiler	GeoMx* Digital Spatial Profiler Analysis Instrument. Includes 1 year manufacturers warranty.	GMX-DSP-1Y
	GeoMx* Digital Spatial Profiler Analysis Instrument. Includes 1 year manufacturers warranty and 1 year service contract.	GMX-DSP-2Y
	GeoMx* Digital Spatial Profiler Analysis Instrument. Includes 1 year manufacturers warranty and 2 year service contract.	GMX-DSP-3Y
	GeoMx* Digital Spatial Profiler Analysis Instrument. Includes 1 year manufacturers warranty and 3 year service contract.	GMX-DSP-4Y
	GeoMx* Digital Spatial Profiler Analysis Instrument. Includes 1 year manufacturers warranty and 4 year service contract.	GMX-DSP-5Y

Materials Included

Item	Quantity
Logitech MK540 Advanced Keyboard, Mouse set	1
LED Monitor 4K 8MP Color 27"	1
Slide Holder Assembly	1
96 Well Plate Adapter	1
Uninterruptible Power Supply	1

Request a quote: <u>GeoMx@nanostring.com</u>

References

- 1 Gupta S, Zugazagoitia J, Martinez-morilla S, Fuhrman K, Rimm DL. Digital quantitative assessment of PD-L1 using digital spatial profiling. Lab Invest. 2020; epub ahead of print.
- 2 Merritt CR, Ong GT, Church SE et al. Multiplex digital spatial profiling of proteins and RNA in fixed tissue. Nat Biotechnol 2020; 38:586-599.

©2020 NanoString Technologies, Inc. All rights reserved. NanoString, NanoString Technologies, nCounter, and the NanoString logo are trademarks or registered trademarks of NanoString Technologies, Inc., in the United States and/or other countries.