# #8042 A complete pipeline for high-plex spatial proteomic profiling and analysis of neural cell phenotypes on a Spatial Molecular Imager and a Spatial **Informatics Platform.**

Shilah A Bonnett\*1, Alyssa Rosenbloom<sup>1</sup>, Mark Conner<sup>1</sup>, Tiên Phan-Everson<sup>1</sup>, Zachary Lewis<sup>1</sup>, Giang Ong<sup>1</sup>, Yan Liang<sup>1</sup>, Emily Brown<sup>1</sup>, Liuliu Pan<sup>1</sup>, Aster Wardhani<sup>1</sup>, Mithra Korukonda<sup>1</sup>, Carl Brown<sup>1</sup>, Iiuliu Pan<sup>1</sup>, Aster Wardhani<sup>1</sup>, Mithra Korukonda<sup>1</sup>, Carl Brown<sup>1</sup>, Liuliu Pan<sup>1</sup>, Aster Wardhani<sup>1</sup>, Mithra Korukonda<sup>1</sup>, Carl Brown<sup>1</sup>, Iiuliu Pan<sup>1</sup>, Iiuliu Pan<sup>1</sup>, Aster Wardhani<sup>1</sup>, Iiuliu Pan<sup>1</sup>, Iiuliu Pan<sup>1</sup>, Iiuliu Pan<sup>1</sup>, Aster Wardhani<sup>1</sup>, Iiuliu Pan<sup>1</sup>, Aster Wardhani<sup>1</sup>, Iiuliu Pan<sup>1</sup>, Iiuliu Vikram Devgan<sup>1</sup>, Michael Rhodes<sup>1</sup>, Gary Geiss<sup>1</sup>, Joseph M Beechem<sup>1</sup>

<sup>1</sup>NanoString<sup>®</sup> Technologies, Seattle, Washington, USA

Abstract

The brain is complex and heterogeneous where cell function and cell-to-cell communication are critical for rapid and accurate performance. The ability to explore protein-driven activities at high resolution within the spatial context of their immediate environment is critical to gain comprehensive pictures of brain development, activity, aging, disease or dysfunction, and inflammatory responses. Many existing approaches for high-plex single-cell spatial proteomics face issues around simplicity, speed, scalability, and big data analysis.

Here, we present an integrated workflow that addresses key concerns around high-plex proteomics. The CosMx<sup>™</sup> Spatial Molecular Imager (SMI) and AtoMx<sup>™</sup> Spatial Informatics Platform comprise an end-to-end workflow that efficiently handles highly multiplex protein analysis at plex sizes exceeding 68 targets. The CosMx protein assays use oligonucleotideconjugated antibodies, that are detected using universal, multi-analyte CosMx readout reagents. The CosMx Mouse Neural Cell Typing and Alzheimer's Pathology panel is optimized to comprehensively profile neural cell lineages across the brain as well as the progression of Alzheimer's disease (AD). Furthermore, 80% of the antibodies making up both panels are crossreactive with human tissue antigens. The AtoMx spatial informatics platform provides full analysis support, including whole-slide image viewer, and methods for performing built-in or fully customizable analyses for cell typing, ligand-receptor analysis, neighborhood analysis and spatial differential expression.

The CosMx protein assay reagents were validated on the FFPE adult mouse brain, mouse embryo, and Alzheimer's positive human brain. We used the CosMx Mouse Neural Cell Typing and Alzheimer's Pathology panel with the CosMx SMI to identify multiple neuronal subtypes, different reactive states of astrocytes and microglia, cell degeneration and proliferation. In a single-cell exploration of mitochondria, we noted distinct patterning of key immune targets based on their immediate microenvironment. Additionally, evaluated the co-expression patterns and activation states of microglia within 200  $\mu$ m from the amyloid plague and tau tangles.

The CosMx SMI is a high-plex spatial multi-omics platform that enables the detection of > 68 proteins at subcellular resolution. In combination with the high-plex CosMx Mouse Neural Cell Typing and Alzheimer's Pathology panel, we present a flexible and scalable informatics platform, a robust solution for comprehensive neural and disease phenotyping that captures the complexity of neuronal and glial cellular activity with full spatial context.



FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. www.nanostring.com | info@nanostring.com © 2023 NanoString Technologies, Inc.



# Differential expression of microglial markers and phosphorylated Taus across an FOV-elaborate signal of individual cells across an FOV



- Biotechnology.
- spatial transcriptomics. BioRxiv.

# Conclusion

1 GFAP GAD67 P2ry12 ChA

- End-to-end workflow from sample preparation to data analysis
- Enables high-plex proteomics up to 68 proteins from a single FFPE slide with spatial context at single-cell resolution

CD68 TMEM119 GFAP DAP12 P2ry12 NeuN

- Protein assay uses validated antibodies conjugated with oligonucleotides,
- detected via microfluidic-based cyclical immunofluorescence imaging process

### Neural Single Cell Spatial Proteomics on Human Alzheimer's Diseased Tissue

# Evaluate the pathology of co-expression patterns of phosphorylated Taus and the activation states of microglia in proximity of amyloid Beta plagues





## References

He et al., 2022. High-plex imaging of RNA and proteins at subcellular resolution in fixed tissue by spatial molecular imaging. Nature

2. Danaher et al., 2022. Insitutype: likelihood-based cell typing for single cell

Uhlén M et al., 2015. Tissue-based map of the human proteome. Science.

Scan here to download or learn more

