




GeoMx[®] IO Proteome Atlas

Spatial Proteomics Redefined

Comprising nearly all Abcam's IHC-validated human antibodies for immuno-oncology (IO) research, the GeoMx IO Proteome Atlas (IPA) is the highest-plex panel available for spatial proteomics of FFPE and fresh frozen tissue sections and enables broad biomarker discovery for translational research. Rapidly discover new protein biomarkers and drug targets by non-destructively profiling over 570 proteins in different tissue compartments such as the tumor, microenvironment, and the immune infiltrate.




Product Highlights

- Find the next blockbuster IO drug target
- Be the first to discover new multiomic biomarkers utilizing spatial proteogenomics with the GeoMx[®] Whole Transcriptome Atlas (WTA)
- Screen the world's largest catalog of IHC-validated antibodies on a single slide
- Customize with up to 40 additional protein targets of your choice

Abcam's IO Catalog at your Fingertips

Find the next blockbuster IO spatial biomarker by profiling the expression of 570+ human protein targets on one slide. Take advantage of the open-source clone information from [Abcam](#) for downstream validation. Although the content is ideal for immuno-oncology research with comprehensive coverage of all 10 Hallmarks of Cancer, the GeoMx IO Proteome Atlas has broad applicability for other research areas such as autoimmunity, infectious disease, and organ transplant.



Content breakdown as distributed across the 10 Hallmarks of Cancer

Flex your Protein Muscles. All 570+ of Them.

Rapidly discover new protein biomarkers and drug targets from human samples with unrivaled plex, getting more comprehensive biological information per slide than traditional IHC or IF assays. Screen multiple samples per slide with whole slide scanning or maximize the number of samples per slide with a Tissue Microarray (TMA). Work with a variety of sample types including whole tissue sections and organoids. Take advantage of NGS readout and use higher capacity Illumina sequencers such as the NovaSeq or NextSeq series with the GeoMx IPA.

Panel Annotations

Content in the GeoMx IPA can be divided into 77 annotated pathways across 556 genes and includes coverage of 60+ phosphorylated targets.

Pathway	# of Targets
AMPK Signaling	33
Androgen Signaling	28
Apoptosis	82
Autophagy	51
Cancer Antigens	2
Cell Adhesion & Motility	102
Cell Cycle	92
Chemokine Signaling	56
Complement System	11
Cytotoxicity	3
Differentiation	125
DNA Damage Repair	36
Drug Resistance	1
EGFR Signaling	8
EMT	95
Endocytosis	38
Epigenetic Modification	82
ERBB2 Signaling	16
Estrogen Signaling	62
FGFR Signaling	13
FoxO Signaling	57

Pathway	# of Targets
Glutamine Metabolism	1
Glycolysis & Glucose Transport	9
Hedgehog Signaling	12
HIF1 Signaling	56
Hippo Signaling	5
IDH1/2	6
IL-1 Signaling	37
IL-17 Signaling	28
IL-2 Signaling	23
IL-6 Signaling	7
Immortality & Stemness	24
Insulin Signaling	43
Interferon Response Genes	10
JAK-STAT Signaling	49
Lipid Metabolism	22
Lymphocyte Regulation	36
Lymphocyte Trafficking	27
Lysosome	9
MAPK Signaling	108
Matrix Remodeling and Metastasis	18
MET Signaling	14

Pathway	# of Targets
MHC Class I Antigen Presentation	24
MHC Class II Antigen Presentation	6
Mitochondrial Metabolism / TCA	15
mTOR Signaling	61
Myc	14
Myeloid Inflammation	43
NF-κB Signaling	45
NO Signaling	7
Notch Signaling	53
NRF2 Signaling	15
Other Interleukin Signaling	112
Oxidative Stress	99
p53 Signaling	58
PDGF Signaling	20
PI3K-Akt Signaling	119
PPAR Signaling	4
Prostaglandin Inflammation	4


Pathway	# of Targets
Proteotoxic Stress	9
Purinergic Signaling	1
RAGE Signaling	9
Senescence	111
T Cell Exhaustion	11
T Cell Checkpoints	18
TCR Signaling	62
TGF-beta Signaling	40
TH1 Differentiation	12
TH17 Differentiation	24
TH2 Differentiation	14
TLR Signaling	74
TNF Signaling	54
Treg Differentiation	13
Type I Interferon Signaling	16
Type II Interferon Signaling	24
VEGF Signaling	40
Wnt Signaling	60

Functional Validation


All antibodies in the GeoMx IPA have been validated with IHC and a host of other orthogonal single plex methods by Abcam, as well as in multiplex by NanoString, to ensure that every antibody has appropriate specificity, sensitivity, and performance. Multiplex validation was completed using FDA-approved TMAs across 77 normal human cores representing 29 tissue types and 54 tumor human cores representing 26 tissue types. All panel antibodies were additionally validated in multiplex across 90+ human cell lines for specificity and sensitivity.

Single plex IHC validation by Abcam


CD3e (Clone ID: EP449E)



Western Blot



Positive IHC staining in normal spleen tissue




Negative IHC staining in normal brain tissue

Multiplex validation by NSTG

Measured on 90 pellet CPA and FDA-approved screening arrays (TMA) with tumor and normal tissues.

Assessed:

- Clear positive signal in positive control cells & tissues
- Clear negative signal in negative control cells & tissues
- Phospho-specific antibodies also validated using established models, including phosphatase inhibitor treatment prior to fixation




GeoMx images of FDA-approved tissue microarrays (TMAs) used for functional validation of the GeoMx IPA


Enable Spatial Proteogenomics

Be the first to discover new biology and novel multiomic biomarkers by utilizing our Spatial Proteogenomics workflow to run the GeoMx IO Proteome Atlas and the GeoMx Human Whole Transcriptome Atlas (WTA) simultaneously on the same tissue section. Updated probe chemistry for GeoMx IPA enables same-slide co-detection of RNA and protein at the same sequencing depth and cost as separate detection of RNA and protein from serial sections.

Performance Remains Consistent from IPA alone to IPA+WTA:



Similar SNR Distribution Pattern



Great Concordance

Ordering Information

GeoMx protein assays arrive ready-to-use and generally ship within 24 hours following purchase.

Item Number	Product Name	Product Description	Quantity
121300160	GeoMx Human IO Proteome Atlas Kit	A 570+ plex protein profiling panel that contains extensive immune, immuno-oncology, metabolism, and additional targets. 4 slide aliquot.	4 slides
121400207	GeoMx Pro Code Pack: Z, Y	NGS readout reagents for GeoMx DSP Protein analysis. Includes Pro Code Z & Y UDI primer plates, and enzyme master mix for each plate.	2 plates, 190 AOI
100473	GeoMx DSP Collection Plate	Barcoded collection plates for use on the GeoMx DSP. Required for AOI tracking.	4 plates, 384 AOI
100474	GeoMx DSP Instrument Buffer Kit	Instrument buffer kit for GeoMx DSP experiments.	24 slides, 12 AOI/slide
121300312	GeoMx Protein Slide Prep Kit	Slide prep kit for GeoMx Protein experiments.	12 slides
121300301	GeoMx Solid Tumor TME Morphology Kit	Morphology kit for visualization of human solid tumors and the tumor microenvironment. For use with protein assays. Includes fluorescent antibodies against Pan-CK, CD45, and a nuclear stain.	12 slides

To view the annotated protein target list for the GeoMx® IO Proteome Atlas visit nanosttring.com/products/geomx-digital-spatial-profiler/io-proteome-atlas

NanoString Technologies, Inc.

530 Fairview Avenue North
Seattle, Washington 98109

T (888) 358-6266
F (206) 378-6288

nanosttring.com
info@nanosttring.com

Sales Contacts

United States us.sales@nanosttring.com
EMEA: europe.sales@nanosttring.com

Asia Pacific & Japan apac.sales@nanosttring.com
Other Regions info@nanosttring.com

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

©2023 NanoString Technologies, Inc. All rights reserved. NanoString, NanoString Technologies, GeoMx, and the NanoString logo are trademarks or registered trademarks of NanoString Technologies, Inc., in the United States and/or other countries.

OCT 2023 MK6080

