

#6

Interrogation of Spatial Transcriptomics on Canine Tumor and Normal Tissue

Sarah Church, Ben Sutton, Erin Piazza, and Kyla Teplitz
NanoString Technologies, Seattle, WA 98109

Abstract

Combination therapy to treat hematological and solid malignancies, including chemotherapy, radiation, targeted and immunotherapy, all hold huge potential for eliciting clinical responses. Informative pre-clinical testing of these approaches can be greatly facilitated using immune competent animals with spontaneous tumors. Most pet dogs are immunologically outbred, immune competent and develop spontaneous tumors such as non-Hodgkin’s lymphoma, glioblastoma, osteosarcoma, urothelial carcinoma and melanoma that share remarkable clinical, biological and genetic features with their human counterparts. As such, pre-clinical testing of therapeutic approaches in dogs with cancer promises to accurately inform human clinical trial design. For this comparative approach to provide maximum information to accelerate human clinical translation of novel combination therapies and identify correlative biomarkers of therapeutic response, it is necessary to develop research tools for deep interrogation of the canine tumor microenvironment (TME). Here, we present spatial transcriptomic analysis of multiple canine tumor and tissue types using the GeoMx® digital spatial profiler (DSP) Canine Cancer Atlas panel.

GeoMx Canine Cancer Atlas

Formalin-fixed paraffin-embedded (FFPE) slides or tissue microarrays were used to profile tumor and normal tissue from canines. Each slide was stained with tissue specific immunofluorescent antibodies, including pan-cytokeratin (PanCK), CD45, Vimentin, IBA1, and CD3. Regions of interest were selected to assess the TME and normal tissue for a comparison as possible. Tissue slides were then run on the DSP instrument using standard DSP methods. We leveraged the Canine Cancer Atlas panel that contains 1,962-canine-specific genes for transcriptomic profiling.




Figure 1. GeoMx Digital Spatial Profiling with the Canine Cancer Atlas enables quantification of 1,962 protein-coding genes in precisely defined regions of interest. The spike-in up to 400 custom RNA targets is also available to add-in probes for non-coding RNA, synthetic DNA, or exogenous genes.

Methods for Transcriptomic Profiling using the GeoMx Canine Cancer Atlas

