

Application of the NanoString nCounter System as an Alternative Method to Investigate Molecular Mechanisms Involved in Host Plant Responses to *Plasmodiophora brassicae*

Background

Clubroot, caused by *Plasmodiophora brassicae*, is a major disease of Brassica crops, including canola, with significant economic impact in Canada. As new pathotypes of *P. brassicae* continue to emerge and overcome host resistance, there is a growing need to better understand host–pathogen interactions at the molecular level. While RNA-seq and qPCR have been widely used for transcriptomic studies and data validation, these methods can be resource-intensive, particularly when investigating defined gene sets across large sample cohorts. The nCounter platform offers a targeted approach that can streamline these efforts.

Research Questions

- Can the nCounter platform reliably validate RNA-seq data from clubroot-infected plant tissues?
- Which genes are consistently involved in clubroot resistance or susceptibility, and could be targeted in future expression panels?
- How does the performance of nCounter compare to RNA-seq and qPCR in terms of time, cost, and accuracy?

Results & Conclusions

- Expression results from the nCounter system showed strong correlations with both RNA-seq and qPCR ($R > 0.90$, $p < 0.01$), demonstrating its reliability and reproducibility in agricultural samples.
- Genes involved in salicylic acid, jasmonic acid, and auxin signaling pathways were differentially expressed and may contribute to disease resistance phenotypes.
- Candidate markers validated by nCounter include WRKY70, PNP-A, DMR6, CHI, and GH3.12, which are implicated in defense regulation and hormone signaling.
- nCounter successfully profiled 29 genes across multiple time points and cultivars, capturing consistent expression trends with fewer processing steps, lower costs, and less technical variability than qPCR or RNA-seq.
- The robustness and efficiency of results generated with the nCounter platform support its adoption as a single end-to-end solution for medium-scale targeted expression studies.

Experimental Setup

Sample Type	Total RNA extracted from pooled rutabaga root tissue
Tissue Type	Root tissue from resistant and susceptible <i>B. napus</i> cultivars
Assay	nCounter Custom CodeSet targeting 29 genes of interest and 4 reference genes
Analyte	RNA
Instrument	nCounter Analysis System




Figure 1: Correlation plots comparing of gene expression measurements between the nCounter platform and RNA-seq at several timepoints, demonstrating high concordance between the methods.

Figure reproduced with permission from Zhou et al. *International Journal of Molecular Sciences*, 23(24), 15581 (2022) under the [Creative Commons license](https://creativecommons.org/licenses/by/4.0/).

Zhou, Q., Galindo-González, L., Hwang, S. F., & Strelkov, S. E. (2022). Application of the NanoString nCounter system as an alternative method to investigate molecular mechanisms involved in host plant responses to *Plasmodiophora brassicae*. *International Journal of Molecular Sciences*, 23(24), 15581. <https://doi.org/10.3390/ijms232415581>

Bruker Spatial Biology | For more information, visit nanosttring.com

Bruker Spatial Biology, Inc.

530 Fairview Avenue North T (888) 358-6266 nanostring.com
Seattle, Washington 98109 F (206) 378-6288 customerservice.bsb@bruker.com

Sales Contacts

United States nasales.bsb@bruker.com Asia Pacific & Japan apacsales.bsb@bruker.com
EMEA: emeasales.bsb@bruker.com Other Regions globalsales.bsb@bruker.com

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

© 2025 Bruker Spatial Biology, Inc. All rights reserved. NanoString, NanoString Technologies, the NanoString logo, nCounter, GeoMx, CosMx, and AtoMx are trademarks or registered trademarks of Bruker Spatial Biology, Inc., in the United States and/or other countries. Any other trademark that appears in this document is the property of its respective owner.