

nCounter<sup>®</sup>

# RUO Lymphoma Subtyping Test

## User Manual

## Changes in this Revision

- Updates to text for clarity and accuracy.
- Updates to list of required materials.
- Updates to Quick Reference Guide.
- Revised Capture Master Mix instructions (page 12, Step 6) to align with best practices.

## Table of Contents

Introduction.....	3
Workflow .....	3
Materials and Equipment .....	4
Materials Supplied by NanoString.....	4
Additional Materials Required .....	4
Equipment.....	4
Sample Preparation Recommendations .....	5
Tissue Processing.....	5
RNA Extraction and Input to RUO LST Assay .....	5
RUO LST Hybridization Protocol.....	6
Quick Reference for LST RUO Hybridization Protocol .....	8
nCounter Panel Plus Products.....	10
Preparing a Merged Reporter Library File (RLF) .....	10
RUO LST Panel Plus Hybridization Protocol.....	11
Quick Reference for LST RUO Panel Plus Hybridization Protocol.....	13
Technical Support.....	15

## Introduction


The RUO Lymphoma Subtyping Test (LST) CodeSet contains probes for the 20 genes within the LST signature: 15 target genes and 5 reference genes.

The signature was originally identified by performing gene expression profiling of Diffuse Large B-Cell Lymphoma (DLBCL) samples using microarrays (Alizadeh et al, 2000). The gene expression profiles indicated two different stages of B-cells differentiation: germinal center (GCB) and activated B cell based (ABC) based on the cell of origin and distinct mechanisms of oncogenesis. This signature was later refined to 20 genes using the FFPE tissues on the nCounter® Analysis System (Scott et al, 2014).

NanoString developed a multiplex gene expression assay on the nCounter Analysis System that can identify the cell of origin (COO) subtypes (ABC-type or GCB-type) of DLBCL. The 20-gene assay was trained using 51 FFPE biopsies. Comparisons were made with COO assignment using the original COO model on matched frozen tissue. A linear predictor score (LPS) calculated as the weighted sum of the 15 normalized target genes compared to pre-defined thresholds determine the DLBCL cell-of-origin subtype. An LST data analysis report is available for purchase separately through NanoString.

These signatures are for research use only (RUO) and not for diagnostic use.

The prepared RNA samples are used as input for nCounter hybridization reactions containing NanoString Reporter and Capture probes. These overnight hybridization reactions enable specific hybridization of Reporter and Capture probes to their target (**Figure 1**). After hybridization, reactions are purified and imaged on an nCounter system.



**Figure 1.** CodeSet chemistry: Capture and Reporter probes bind to the target.

All nCounter gene expression panels are sold in increments of 12 reactions. Additional materials to process samples on an nCounter Analysis System or SPRINT Profiler are required and sold separately (see [Table 1](#)).

## Workflow

1. Prepare FFPE RNA samples and assess quantity and quality (refer to [Preparing RNA from FFPE Samples \(MAN-10050\)](#)).
2. Set up RNA hybridization to the RUO LST CodeSet (0.5 hours).
3. Hybridize overnight in thermal cycler with a heated lid (15–24 hours).
4. Purify and bind hybridized RNA targets on an nCounter Analysis System (Pro or MAX/FLEX) or SPRINT Profiler.

## Materials and Equipment

### Materials Supplied by NanoString

**Table 1.** NanoString-provided materials required to run nCounter RUO LST Assay.

Item	Reagents	Storage
RUO LST CodeSet - <i>Catalog # RUO-LST-12</i>	RUO LST Reporter CodeSet tube RUO LST Capture ProbeSet tube	-80°C -80°C
RUO LST Reference Sample (optional) - <i>Catalog # NALG-REF-LST</i>	RUO LST Reference Sample	-80°C
RUO LST Panel Plus (optional) - <i>Catalog # varies</i>	RUO LST Reporter Plus CodeSet tube RUO LST Capture Plus CodeSet tube	-80°C -80°C
nCounter Master Kit ( <i>for MAX/FLEX/Pro</i> ) - <i>Catalog # NAA-AKIT-012</i>	nCounter Sample Cartridge Prep Plate Prep Pack, including Hybridization Buffer	-20°C 4°C 15-25°C
nCounter SPRINT Reagent Pack ( <i>for SPRINT</i> ) - <i>Catalog # SPRINT-REAG-KIT</i>	nCounter SPRINT Reagent C nCounter SPRINT Reagent A, B, and Hybridization Buffer	4°C 15-25°C
nCounter SPRINT Cartridge ( <i>for SPRINT</i> ) - <i>Catalog # SPRINT-CAR-1.0</i>	nCounter SPRINT Cartridge	-20°C

### Additional Materials Required

**Table 2.** Additional materials required (not provided by NanoString).

Item	Supplier
Multi-channel pipette for 200 µL	Various
Single-channel pipettes for 2 µL, 20 µL, 200 µL, 1000 µL	
RNase-free pipette tips with aerosol barriers	
0.2-mL strip tubes and caps, nuclease-free (SPRINT users only; these are provided in Master Kits for MAX/FLEX/Pro users)	
Disposable gloves	
Molecular biology-grade nuclease-free water	

### Equipment

**Table 3.** Required equipment to run the nCounter RUO LST Assay.

Equipment	Supplier
Thermal cycler with a programmable heated lid	Various
Standard benchtop centrifuge with a fixed-angle rotor that fits 1.5 mL tubes	
Picofuge or mini-centrifuge with 1.5 mL-tube rotor and strip tube rotor	
NanoString nCounter® Pro, MAX, or FLEX Analysis System or SPRINT Profiler	NanoString

## Sample Preparation Recommendations

### Tissue Processing

NanoString recommends submitting a pathology-reviewed H&E-stained slide and corresponding unstained slides from an FFPE block (meeting the requirements defined in [Table 4](#)) to the lab for testing. A board-certified pathologist should circle the area of viable tumor on the H&E-stained slide to enable macrodissection (if required) and an assessment of tumor cellularity. It is recommended that the tissue sections submitted for testing by LST have at least 10% tumor cellularity of the sample. Macrodissection is desired, when possible, if the tumor constitutes less than 50% of the total sample area, to avoid potential bias towards GCB.

[Table 4](#) lists the number of unstained slide-mounted tissue sections recommended for one RNA extraction based on the tumor surface area ( $\text{mm}^2$ ) estimated by a pathologist using an H&E-stained slide.

**Table 4.** Recommended numbers of unstained slide-mounted tissue sections.

Tumor Surface Area	Number of Slides (5 $\mu\text{m}$ thickness)*	
	Minimum Number of Slides	Recommended Input
2 $\text{mm}^2$	4	8 or more, as available
3 $\text{mm}^2$	3	8
4–7 $\text{mm}^2$	2	5–7
8–15 $\text{mm}^2$	1	3–4
$\geq 16 \text{ mm}^2$	1	2

\*Slide-mounted tissue sections between 4–10  $\mu\text{m}$  thickness are also acceptable. The number of slides should be adjusted to obtain a minimum total tissue volume of 0.04  $\text{mm}^3$ .

### RNA Extraction and Input to RUO LST Assay

See [Preparing RNA from FFPE Samples \(MAN-10050\)](#) for guidance to extract RNA and determine its concentration and purity.

Isolated RNA is expected to meet the following specifications:

- RNA concentration:  $\geq 6.25 \text{ ng}/\mu\text{L}$
- RNA purity: A260/A280 ratio between 1.65 and 2.35

The recommended RNA input for this assay is 500 ng. The acceptable RNA input range for hybridization is 62.5–1000 ng.

- Calculate the volume (in microliters) of RNA sample to add to the hybridization reaction by dividing the desired sample input (e.g., 500 ng) by the measured concentration.
- If the calculated concentration of the sample is between 6.25  $\text{ng}/\mu\text{L}$  and 50  $\text{ng}/\mu\text{L}$ , add the maximum volume of 10  $\mu\text{L}$ .
- For samples that require less than 10  $\mu\text{L}$ , calculate the volume of water required to generate a 10  $\mu\text{L}$  final sample volume.

## RUO LST Hybridization Protocol

### Reporter CodeSet and Capture ProbeSet Handling Instructions:

- During setup, do not vortex or pipette vigorously to mix. Instead, gently flick or invert the tubes.
- To spin down contents of tubes, a picofuge or mini-centrifuge is recommended. If using a centrifuge, spin at <3000xg for <10 seconds. Do not “pulse” spin as it will cause the centrifuge to go to maximum speed and may spin the probes out of solution.

**IMPORTANT:** Check the reagent labels before you begin to ensure use of the correct reagents. If using RUO LST Panel Plus, refer to [RUO LST Panel Plus Hybridization Protocol](#) on page 11.

1. **Pre-heat** a thermal cycler to **65°C** with a heated lid at **70°C**; set the calculated reaction volume to 30  $\mu\text{L}$  and the time interval to “infinite”. **Do not** set the thermal cycler to ramp down to 4°C at the end of the run.

**NOTE:** A thermal cycler with a heated lid is required for this protocol. NanoString recommends a thermal cycler with a *programmable* heated lid. Models without programmable lids, if used, should be set to ensure that the heated lid does not exceed 110°C.

2. **Remove** Reporter CodeSet and Capture ProbeSet reagents from the -80° freezer and thaw at room temperature, shielded from light. Once thawed, invert or gently flick the tubes several times to mix well, then briefly spin down reagents. Store reagents on ice if not proceeding immediately.

**IMPORTANT:** After it has thawed, inspect the tube of Reporter CodeSet to make sure no colored precipitate is present. If you see a colored precipitate, heat the entire tube to 75°C for 10 minutes and cool at room temperature before using.

3. **Create a hybridization master mix** by adding the Hybridization Buffer to the tube containing the Reporter CodeSet at room temperature. Do not remove the Reporter CodeSet from the tube; add the Hybridization Buffer directly into the CodeSet tube. **Do not add the Capture ProbeSet to the master mix.**

**Table 3.** LST RUO Master Mix for one nCounter assay (12 reactions + 1 reaction dead volume).

Component	CodeSet Master Mix ( $\mu\text{L}$ )	Per Reaction ( $\mu\text{L}$ )
Reporter CodeSet	65 (in tube)	5
Hybridization Buffer	130	10
<b>Total Volume</b>	<b>195</b>	<b>15</b>

4. **Gently flick or invert the hybridization master mix tube repeatedly** to mix, then briefly spin down.
5. **Label a strip tube.** If necessary, cut strip in half to fit in a picofuge with strip tube adaptor, and label both halves. For MAX/FLEX/Pro users, use the strip tubes provided with the nCounter Master Kits, ensuring that the notch is positioned between tubes 1-2 and 8-9.


6. **Prepare hybridization reactions** using a new pipette tip at every step:
  - a. **Add 15  $\mu\text{L}$  of hybridization master mix** to each well of the prepared strip tube.
  - b. **Add up to 10  $\mu\text{L}$  of sample** to each tube containing hybridization master mix. If using less than 10  $\mu\text{L}$  of sample, add RNase-free water to each tube to bring the volume to 25  $\mu\text{L}$ .
  - c. **Mix the Capture ProbeSet tube** by inverting or flicking, and briefly spin down the contents.
  - d. **Add 5  $\mu\text{L}$  of Capture ProbeSet** to each tube.
  - e. **Cap the strip tubes tightly and mix** by inverting the tubes several times and gently flicking to ensure complete mixing.
  - f. **Spin briefly and immediately place** the tubes in a pre-heated **65°C** thermal cycler.
7. **Incubate hybridization reactions at 65°C for 15–24 hours.** Hybridizations should be left at 65°C until ready for processing.
8. Once the hybridization reactions have been removed from the thermal cycler, **proceed immediately to an nCounter Prep Station or SPRINT Profiler.** Please refer to instrument-specific user manuals ([nCounter Pro Analysis System User Manual \(MAN-10147\)](#), [nCounter Analysis System User Manual for MAX/FLEX Systems \(MAN-C0035\)](#), [nCounter SPRINT Profiler User Manual \(MAN-10017\)](#)).

# RUO LST Hybridization Protocol

Please read through the full protocol before beginning. This illustrated workflow is intended for quick reference at the bench.


**1 Prepare for hybridization**

- Preheat thermocycler to **65° C** with a heated lid at **70° C**.
- Thaw** Codeset and samples.



**2 Create & aliquot Hybridization Master Mix**


- Add **130 µL of Hybridization Buffer to the Reporter CodeSet tube** to create Hybridization Master Mix.
- Flick to mix, then briefly spin down contents.
- Aliquot **15 µL of Hybridization Master Mix** into each tube of a labeled 12-tube strip.



**3 Add sample**


- Add **up to 10 µL of sample** to each tube.

*Note: If using less than 10 µL of sample, add RNase-free water to each tube to bring the volume to 25 µL*



**4 Add Capture ProbeSet**


- Flick-mix** Capture ProbeSet and **spin down briefly**.
- Add **5 µL of Capture ProbeSet** to each tube.
- Cap tightly, **flick-mix**, and **spin down briefly**.





## 5 Hybridize

- ❑ **Incubate** hybridization reactions at **65° C for 15-24 hours**. Hybridizations should be left at 65° C until ready for processing.



### OPTION A: nCounter® Pro or MAX/FLEX



## 6a Load nCounter Pro or MAX/FLEX

See **MAN-C0035, nCounter Analysis System User Manual** or **MAN-10147, nCounter Pro Analysis System User Manual** for details.

### OPTION B: nCounter® SPRINT



## 6b Load nCounter SPRINT

See **MAN-10017, nCounter SPRINT Profiler User Manual** for details.

## nCounter Panel Plus Products

All off-the-shelf nCounter Gene Expression panels are customizable by adding 6 to 55 user-defined probes. These additional probes are referred to as a Panel Plus product. Similarly, customized CodeSets can be supplemented with CodeSet Plus products.

### Preparing a Merged Reporter Library File (RLF)

All nCounter Plus products are accompanied by an add-in library file (ALF) that specifies the association between each Plus probe pair and its target. Information from the ALF must be merged with the reporter library file (RLF) from the CodeSet the Plus product is being added to, prior to scanning the cartridge or the barcodes on the nCounter Digital Analyzer or SPRINT Profiler. Failure to merge an ALF with an nCounter CodeSet RLF will result in no count information being collected for targets of Plus products.

To obtain a merged RLF file, email NanoString at [bioinformatics@nanosttring.com](mailto:bioinformatics@nanosttring.com). Include both the ALF for your Plus product and the RLF for the CodeSet into which you will spike the Plus product. A new merged RLF will be generated and emailed to the requestor that contains all probe information for both the Plus product and the original CodeSet.

**IMPORTANT:** When using a Plus product, you **MUST** use a merged RLF to ensure counting of Plus targets. Ensure that you have the merged RLF file before beginning the hybridization protocol.

## RUO LST Panel Plus Hybridization Protocol

### Reporter CodeSet and Capture ProbeSet Handling Instructions:

- During setup, do not vortex or pipette vigorously to mix. Instead, gently flick or invert the tubes.
- To spin down contents of tubes, a picofuge or mini-centrifuge is recommended. If using a centrifuge, spin at <3000xg for <10 seconds. Do not “pulse” spin as it will cause the centrifuge to go to maximum speed and may spin the probes out of solution.

**IMPORTANT:** Check the reagent labels before you begin to ensure use of the correct reagents. If you are NOT using Panel Plus, refer to the [RUO LST Hybridization Protocol](#) on page 6.

**IMPORTANT:** When using a Panel Plus product, you MUST use a merged RLF to ensure counting of Panel Plus targets. Ensure that you have the merged RLF file before beginning the hybridization protocol. See [Preparing the Merged RLF](#) on page 10 for more information.

1. **Pre-heat** a thermal cycler to **65°C** with a heated lid at **70°C**, set the calculated reaction volume to 30  $\mu\text{L}$ , and the time interval to “infinite”. **Do not** set the thermal cycler to ramp down to 4°C at the end of the run.

**NOTE:** A thermal cycler with a heated lid is required for this protocol. NanoString recommends a thermal cycler with a *programmable* heated lid. Models without programmable lids, if used, should be set to ensure that the heated lid does not exceed 110°C.

2. **Remove** Reporter CodeSet, Capture ProbeSet, Reporter Plus, and Capture Plus tubes from the -80° freezer and thaw at room temperature, shielded from light. Once thawed, invert or flick the tube several times to mix well, then briefly spin down reagents. Store reagents on ice if not proceeding immediately.

**IMPORTANT:** After it has thawed, inspect the tube of Reporter CodeSet to make sure no colored precipitate is present. If you see a colored precipitate, heat the entire tube to 75°C for 10 minutes and cool at room temperature before using.

3. **Create a hybridization master mix** by adding the Reporter Plus CodeSet and Hybridization Buffer to the tube containing the Reporter CodeSet at room temperature. Do not remove the Reporter CodeSet from the tube - add components directly into the CodeSet tube. **Do not add the Capture ProbeSet to the hybridization master mix.**

**Table 4.** RUO LST Hybridization Master Mix for one nCounter assay (12 reactions +1 reaction dead volume).

Component	Hybridization Master Mix ( $\mu\text{L}$ )	Per Reaction ( $\mu\text{L}$ )
Reporter CodeSet	65 (in tube)	5
Reporter Plus CodeSet	26	2
Hybridization Buffer	130	10
<b>Total Volume</b>	<b>221</b>	<b>17</b>

4. **Flick or invert the hybridization master mix tube repeatedly** to mix, **then briefly spin down**.
5. **Label** a strip tube. If necessary, cut strip in half to fit in a picofuge with strip tube adaptor, and label both halves. For MAX/FLEX/Pro users, use the strip tubes provided with the nCounter Master Kits, ensuring that the notch is positioned between tubes 1-2 and 8-9.
6. **Prepare hybridization reactions** using a new pipette tip at every step:
  - a. **Add 17  $\mu$ L of hybridization master mix** to each well of the prepared strip tube.
  - b. **Add up to 10  $\mu$ L of sample** to each tube containing hybridization master mix. If using less than 10  $\mu$ L of sample, add RNase-free water to each tube to bring the volume to 27  $\mu$ L.
  - c. **Mix the Capture ProbeSet and Capture Plus tubes** by inverting or flicking, and briefly spin down the contents.
  - d. **Create a Capture Master Mix** by adding **14  $\mu$ L** of the Capture Plus to the Core Capture ProbeSet tube. Mix by inverting or gently flicking, and briefly spin down the contents.
  - e. **Add 6  $\mu$ L** of Capture Master Mix to each tube.  
**NOTE:** Final hybridization volume for **RUO LST CodeSet + Panel Plus** is **33  $\mu$ L**.
  - f. **Cap the strip tubes tightly and mix** by inverting the tubes several times and gently flicking to ensure complete mixing.
  - g. **Spin briefly and immediately place** the tubes in a pre-heated **65°C** thermal cycler.
7. **Incubate hybridization reactions at 65°C for 15–24 hours**. Hybridizations should be left at **65°C** until ready for processing.
8. Once the hybridization reactions have been removed from the thermal cycler, **proceed immediately to an nCounter Prep Station or SPRINT Profiler**. Please refer to instrument-specific user manuals ([nCounter Pro Analysis System User Manual \(MAN-10147\)](#), [nCounter Analysis System User Manual for MAX/FLEX systems \(MAN-C0035\)](#), [nCounter SPRINT Profiler User Manual \(MAN-10017\)](#)).

# RUO LST Panel Plus Hybridization Protocol

Please read through the full protocol before beginning. This illustrated workflow is intended for quick reference at the bench.

## 1 Prepare for hybridization


- Preheat thermocycler to **65° C** with a heated lid at **70° C**.
- Thaw** Codeset and CodeSet Plus tubes and samples.



Preheat to 65° C with a heated lid at 70° C

## 2 Create & aliquot Hybridization Master Mix


- Add **130 µL of Hybridization Buffer** and **26 µL of Reporter Plus** to the **Reporter CodeSet** tube to create **Hybridization Master Mix**.
- Flick to mix, then briefly spin down contents.
- Aliquot **17 µL of Hybridization Master Mix** into each tube of a labeled 12-tube strip.



## 3 Add sample


- Add **up to 10 µL of sample** to each tube.

*Note: If using less than 10 µL of sample, add RNase-free water to each tube to bring the volume to 27 µL.*



## 4 Add Capture Master Mix

- Flick-mix** Capture ProbeSet and Capture Plus tubes and **spin down briefly**.
- Add **14 µL of Capture Plus** to the **Capture ProbeSet** to create the **Capture Master Mix**. Flick to mix, then **briefly spin down** contents.
- Add **6 µL of Capture Master Mix** to each tube.
- Cap tightly, **flick-mix**, and **spin down briefly**.




**5 Hybridize**

**Incubate** hybridization reactions at **65° C for 15-24 hours.** Hybridizations should be left at 65° C until ready for processing.



**OPTION A: nCounter® Pro or MAX/FLEX**

**OPTION B: nCounter® SPRINT**



**6a Load nCounter Pro or MAX/FLEX**

See **MAN-C0035, nCounter Analysis System User Manual** or **MAN-10147, nCounter Pro Analysis System User Manual** for details.

**6b Load nCounter SPRINT**

See **MAN-10017, nCounter SPRINT Profiler User Manual** for details.

## Technical Support

For technical support, please contact [support@nanosttring.com](mailto:support@nanosttring.com).

### References Cited

Alizadeh, AA. (2000). Distinct types of diffuse large B-cell lymphoma identified by gene expression profiling. *Nature*, 403, 503-511. <https://www.nature.com/articles/35000501>

Scott, DW. (2014). Determining cell-of-origin subtypes of diffuse large B-cell lymphoma using gene expression in formalin-fixed paraffin-embedded tissue. *Blood*, 123(8), 1214-1217. <https://ashpublications.org/blood/article/123/8/1214/32828/Determining-cell-of-origin-subtypes-of-diffuse>

### Intellectual Property Rights

This nCounter RUO Lymphoma Subtyping Test Assay User Manual and its contents are the property of NanoString Technologies, Inc. ("NanoString"), and are intended for the use of NanoString customers solely in connection with their operation of an nCounter Analysis System or SPRINT Profiler. The nCounter Analysis System and SPRINT Profiler (including both software and hardware components) and this User Manual and any other documentation provided to you by NanoString in connection therewith are subject to patents, copyright, trade secret rights, and other intellectual property rights owned by or licensed to NanoString. No part of the software or hardware may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into other languages without the prior written consent of NanoString. For a list of applicable patents, see [www.nanosttring.com/company/patents](http://www.nanosttring.com/company/patents).

### Trademarks

NanoString, NanoString Technologies, the NanoString logo, nCounter, and nSolver are trademarks or registered trademarks of NanoString Technologies, Inc., in the United States and/or other countries. All other trademarks and/or service marks not owned by NanoString that appear in this document are the property of their respective owners.

© 2022 NanoString Technologies, Inc. All rights reserved.

For Research Use Only. Not for use in diagnostic procedures.

---

**For more information, please visit [nanosttring.com](https://nanosttring.com)**

#### NanoString Technologies, Inc.

530 Fairview Avenue North  
Seattle, Washington 98109

T (888) 358-6266  
F (206) 378-6288

[info@nanosttring.com](mailto:info@nanosttring.com)

#### Sales Contacts

United States: [us.sales@nanosttring.com](mailto:us.sales@nanosttring.com)  
EMEA: [europa.sales@nanosttring.com](mailto:europa.sales@nanosttring.com)

Asia Pacific & Japan: [apac.sales@nanosttring.com](mailto:apac.sales@nanosttring.com)  
Other Regions: [info@nanosttring.com](mailto:info@nanosttring.com)