



## Quick Start Guide

# nCounter® Vantage 3D™ RNA:Protein Solid Tumor Assay for lysate


This quick start guide provides an overview of the research protocol for preparation of RNA and protein from cell or tissue lysate samples. If you are a first-time user and for buffer preparation, please read the full protocols related to your Vantage 3D Assay and use this as a reference in subsequent experiments. Contact NanoString Support ([support@nanostring.com](mailto:support@nanostring.com)) to receive additional assistance with this assay.

### 1 Prepare Sample Lysate (See MAN-10054)

1. Lyse sample by adding appropriate volume of SDS lysis buffer to achieve approximately 0.5 - 1.5 mg/mL protein concentration
2. Boil SDS lysate at **95°C** for **10 minutes**
3. Cool SDS lysate on ice, spin briefly, vortex, and spin briefly again
4. Remove detergent using the Pierce Detergent Removal Spin Columns




5. Determine protein concentration using the Pierce 660nm Protein Assay



6. After determining the protein concentration, samples can be stored at **-80°C** if desired

## 2 Prepare RNA Lysate (See MAN-10054)


1. Dilute a small aliquot of the detergent-free lysate to 0.25 mg/mL using nuclease-free dH<sub>2</sub>O. (A volume of 4 µl is required for hybridization)
2. Transfer diluted lysate to a 12-well strip tube



3. Store at **-80°C** until ready to proceed with hybridization

## 3 Bind Protein to Plate (See MAN-10054)


1. Dilute the detergent-free lysate to 5 µg/ml using PBS pH 7.2 with 2 mM NaN<sub>3</sub>
2. Add 50 µl of the diluted lysate into a well of a protein-binding plate (MAXISORP plate) and seal with parafilm
3. Incubate for at least **2 hours** or **overnight** at room temperature




4. Add 200 µl of Blocking Buffer that has been pre-warmed to room temperature
5. Incubate for **5 minutes** at room temperature (Maximum incubation time is **1 hour**)
6. Remove and discard supernatant by flicking the plate into a sink and striking on a fresh paper towel hard enough to remove any residual liquid
7. Wash 3X by adding 250 µl of TBST pre-warmed to room temperature and incubating for **1 minute**. Remove and discard supernatant by flicking the plate into a sink and striking on a fresh paper towel hard enough to remove any residual liquid

#### 4 Prepare Protein Lysate (See MAN-10054)


1. Prepare a working antibody solution by adding 16  $\mu$ l of the NanoString antibody mix to 625  $\mu$ l of blocking buffer
2. Add 50  $\mu$ l of the working antibody solution to each well
3. Seal the plate and incubate at room temperature for **2 hours** on an orbital shaker at 350-400 RPM (Maximum incubation time is **3 hours**)



4. Carefully remove and discard all supernatant with a single-channel pipette
5. Wash 6X by adding 250  $\mu$ l of TBST pre-warmed to room temperature and incubating for **5 minutes**. Remove and discard supernatant by flicking the plate into a sink and striking on a fresh paper towel hard enough to remove any residual liquid
6. Add 50  $\mu$ l of buffer RLT to each well
7. Seal the plate and incubate at room temperature for **5 minutes** on an orbital shaker at 350-400 RPM



8. Transfer the lysate to a 12-well strip tube



9. See appropriate manual for hybridization set up:  
 MAN-10059 Protein Only  
 MAN-10060 RNA:Protein (CodeSet)  
 MAN-10065 RNA:Protein (TagSet)



---

For more information, please visit [nanosttring.com](https://nanosttring.com)

**NanoString Technologies, Inc.**

530 Fairview Avenue North  
Seattle, Washington 98109

T (888) 358-6266  
F (206) 378-6288

[nanosttring.com](https://nanosttring.com)  
[support@nanosttring.com](mailto:support@nanosttring.com)

**Additional Customer Resources**

For more comprehensive information, visit us at [nanosttring.com](https://nanosttring.com) and go to **Support > Product Support** to view manuals and other technical product literature.

**FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.**

©2017-2019 NanoString Technologies, Inc. All rights reserved. NanoString, NanoString Technologies, the NanoString logo, nCounter and Vantage 3D are trademarks or registered trademarks of NanoString Technologies, Inc., in the United States and/or other countries.