

# MicroRNA Biomarker Discovery and Validation

Using the nCounter® Platform

MK1187 | February 2019 NanoString Technologies®, Inc.

Authors: Lucas Dennis, Michael Rhodes, Kirsteen Maclean, and Anisha Kharkia



# Introduction

MicroRNAs (miRNAs) are a class of small, non-coding RNA that regulate the gene expression of target mRNAs via posttranscriptional gene silencing<sup>1</sup>. These short RNAs have been implicated in the widespread control of critical biological processes such as proliferation, differentiation, and apoptosis<sup>2-4</sup>. Due to their central role in developmental processes, perturbations in miRNA expression patterns have been implicated in many diseases, including cancer, Alzheimer's disease, heart disease, and diabetes<sup>5-7</sup>. Recent work has focused on the promise of developing miRNA expression signatures as prognostic biomarkers of disease <sup>8-11</sup>.

The nCounter miRNA expression essay enables users to rapidly and efficiently profile the top 800 highly-curated human miRNAs from miRBase 22 in a single tube. This product provides a sensitive, reproducible, and highly-multiplexed method for detecting specific miRNAs within purified RNA isolated from any source, including biological fluids, formalin-fixed paraffin-embedded (FFPE) samples, and fresh or frozen cells and tissues<sup>12</sup> (FIGURE 1). Refer to the Tech Note for nCounter miRNA Expression Analysis in Plasma and Serum Samples for additional details. A significant advantage of the assay is that it uses NanoString's nCounter platform to offer direct digital counts of each miRNA without the use of reverse transcription or amplification, making it highly specific. In addition to high performance on FFPE and fresh/frozen tissue, the nCounter miRNA expression assay is highly reproducible (FIGURE 2) with replicate count correlations ( $\mathbb{R}^2$ ) greater than 0.99. It is a unique biomarker discovery and signature development tool that enables collection of expression data (FIGURE 3) in a short amount of time with minimal hands-on manipulation.

nCounter technology is based on a novel method of direct molecular barcoding and digital detection of target molecules through the use of color-coded probe pairs<sup>13</sup>. The nCounter miRNA Sample Preparation Kit provides reagents for ligating unique oligonucleotide tags (miRtags) onto the 3' end of target miRNAs such that short RNA targets can be detected by nCounter probes (FIGURE 4). Sample preparation involves multiplexed ligation of the specific tags to their target miRNA and an enzymatic purification to remove nonligated tags. nCounter miRNA expression assays are highly specific, and can distinguish between miRNAs that differ by just a few base pairs, as shown in cross-hybridization data with the let-7 family, a family of highly homologous miRNAs which can differ in sequence by as little as a single nucelotide (FIGURE 5). Sequence specificity of the ligation reaction is ensured by the use of T<sub>m</sub>-optimized bridging oligos that are complementary to a portion of both the target miRNA and miRNA-specific tag along with careful, stepwise control of annealing and ligation temperatures.



**FIGURE 1:** Expression correlation is preserved when profiling matched FFPE and fresh/frozen tissues. Sections of tissue from a single human liver were formalin-fixed and paraffin-embedded (FFPE) or frozen. 100 ng of total RNA was isolated from the matched FFPE and frozen tissue and profiled using the nCounter miRNA Expression Assay. Count correlations between the samples were high ( $R^2 > 0.95$ ) indicating that generation and handling of FFPE samples has no significant impact on performance of the nCounter assay.



**FIGURE 2:** nCounter Human v3 miRNA Expression assay is highly reproducible. Count correlations between the preparation of the same RNA sample (100 ng) were extremely high (R2 > 0.99) demonstrating the precision of the nCounter Human v3 miRNA assay.



**FIGURE 3:** Differential miRNA expression analysis on the nCounter platform. nCounter miRNA assays are capable of accurately quantifying changes in miRNA expression. Human Reference RNA (100 ng per assay) and Brain Reference RNA (100 ng per assay) were profiled in triplicate using the nCounter Human v3 miRNA assay. Hundreds of miRNAs are significantly differentially expressed between samples.



This approach specifically captures all mature miRNAs in a sample despite the large overall sequence diversity and associated broad  $\rm T_m$  range between miRNAs

(FIGURE 6). Control RNA included in the Sample Preparation Kit is used to monitor ligation efficiency and specificity through each step of the reaction (TABLE 1). All nCounter assays are run using an nCounter instrument. miRNA expression data can be easily imported into the free nSolver<sup>™</sup> Analysis Software

(www.nanostring.com/products/analysis-software/nsolver) for data normalization and visualization, including heat maps, histograms, and violin plots. The nCounter platform has been used to profile miRNAs from a variety of sample types and in a wide range of biological areas of research, including biomarker discovery<sup>14</sup>, signature development<sup>15</sup>, immuno-oncology<sup>16</sup>, oncology<sup>17</sup>, neuroscience<sup>18</sup>, stem cell research<sup>19</sup>, autophagy<sup>20</sup>, infectious disease<sup>21</sup>, immunology<sup>22</sup>, and many others<sup>23,24</sup>.

# **Designed for Biomarker Discovery and Validation**

nCounter miRNA assay content is based on miRBase, a bioinformatics repository for small RNA sequence and annotation information<sup>25,26</sup>. These assays are periodically revised to account for updates to miRBase and to ensure that they include the most biologically-relevant miRNAs and controls. The miRNAs included in the nCounter Human v3 miRNA assay account for 98% of all observed sequencing reads in miRBase 22 (released on 12 March 2018). Additionally, NanoString uses proprietary metrics such as observed read ratios and expression analytics in order to screen potential content for inclusion in the assay. Together, these methods help ensure that the miRNA assay content is heavily weighted toward biologically-significant miRNAs. In addition to this data-driven selection, NanoString also performs literature reviews to ensure that actionable and clinically relevant miRNAs are included in our nCounter miRNA assays<sup>28,29</sup>. This metric is important as not all potentially clinically relevant miRNA biomarkers are classified as "high confidence" by miRBase. NanoString's assay design philosophy uses a holistic set of selection criteria to

Α							
miRNA	Sequence						
Let-7a	UGAGGUAGUAGGUUGUAUAGUU						
Let-7b	UGAGGUAGUAGGUUGU <mark>G</mark> UGGUU						
Let-7c	UGAGGUAGUAGGUUGUAU <mark>G</mark> GUU						
Let-7d	AGAGGUAGUAGGUUG <b>C</b> AUAGUU						
Let-7e	UGAGGUAG <mark>G</mark> AGGUUGUAUAGUU						
Let-7f	UGAGGUAGUAG <mark>A</mark> UUGUAUAGUU						
Let-7g	UGAGGUAGUAG <mark>U</mark> UUGUAUAGUU						
Let-7i	UGAGGUAGUAG <mark>G</mark> UUGU <mark>GCU</mark> GUU						

ensure that each miRNA assay contains a comprehensive collection of miRNAs enriched for biological activity and ideal for both targeted discovery and validation experiments. For greater detail about the Human v3 miRNA expression assay content, refer to the target list (gene list) available on NanoString's website at www.NanoString.com/products/miRNA

A large number of control probes are included in each nCounter miRNA assay to help ensure robust performance. miRNA assay controls have been updated to include probes for both hybridization and ligation based on sequences generated by the External RNA Controls Consortium (ERCC), an independent consortium of academic, public, and private organizations sponsored by the National Institute of Standards and Technology (NIST)<sup>30</sup>. ERCC sequences were developed to enable generation of reproducible research control probes for gene expression analysis and related biological questions via the use of synthetically-derived sequences that do not occur in most genomes. **TABLE 1** outlines each type of control included in nCounter miRNA expression assays.



**FIGURE 4:** Specific capture of miRNA targets via ligation. miRNAs (red) are specifically ligated (circled region) to unique tags (blue) for downstream detection via hybridization with the nCounter miRNA expression assay. Ligation specificity is ensured by  $T_m$ -optimization and  $T_m$ - balancing of the bridge (yellow) sequence that is complementary to each mature miRNA and cognate miRtag coupled with careful thermal control of the ligation reaction.

	Let-7a	Let-7b	Let-7c	Let-7d	Let-7e	Let-7f	Let-7g	Let-7i
Let-7a	100%	1%	5%	4%	17%	4%	0%	0%
Let-7b	0%	100%	0%	0%	0%	0%	0%	0%
Let-7c	0%	11%	100%	0%	0%	0%	0%	0%
Let-7d	0%	0%	0%	100%	0%	0%	0%	0%
Let-7e	0%	0%	0%	0%	100%	0%	0%	0%
Let-7f	1%	0%	0%	0%	0%	100%	0%	0%
Let-7g	0%	0%	0%	0%	0%	0%	100%	0%
Let-7i	0%	0%	0%	0%	0%	0%	0%	100%

**FIGURE 5:** nCounter miRNA Expression assays are highly specific. (A) Single-nucleotide mismatch cross-hybridization rates were empirically determined for the let-7 family, a family of highly homologous miRNAs which can differ in sequence by as little as a single nucleotide. (B) Each let-7 miRNA target was assayed alone and in the presence of all let-7 probes. Cross-hybridization was calculated as: (off-target counts/on-target counts).





FIGURE 6: Calculated Tm distribution of human mature miRNA sequences. Percentage G/C content of mature miRNAs is highly variable, causing a broad Tm distribution. The nCounter miRNA sample preparation method normalizes this Tm distribution while preserving miRNA specificity via miRtag ligation in order to generate ideal targets for downstream detection via hybridization with nCounter Capture and Reporter probes.

# miRNA assay Performance and Application Data

One of the benefits of the direct digital counting methodology afforded by the NanoString nCounter platform is that it is perfectly positioned for translational research needs, enabling the utility of very poor quality RNA, particularly from highly-degraded clinical FFPE samples to address novel treatment options in cancer<sup>31</sup>. Profiling of miRNAs from FFPE and their potential as novel biomarkers in a variety of tumor types has therefore gained much traction. Recent examples have focused on the identification of novel miRNAs in the progression of colorectal cancer<sup>32,33</sup>. Using both fresh frozen and FFPE human tissue as well as mouse models of colorectal cancer, Valeri et al. performed nCounter miRNA profiling on extracted RNA and characterized the up-regulation of miR-135b as a potential biomarker for colorectal cancer<sup>32</sup>. Furthermore, a comprehensive study by Hur et al. identified a metastatic-specific miRNA signature in patient colorectal cancers, providing evidence that miRNAs may be clinically applicable to predict prognosis and distant metastasis in colorectal cancer<sup>33</sup>. To date, triple-negative breast cancer (TNBC) has been one of the most challenging breast tumor subtypes to treat, not least because only a small proportion of patients respond to neoadjuvant therapy regimes<sup>34</sup>. miRNAs have therefore emerged as attractive candidates for novel molecular biomarkers in the progression of breast cancer. Cascione et al. wished to stratify TNBC using NanoString miRNA expression profiling by investigating the profiles in tumor, adjacent, non-tumor,

#### TABLE 1.

Description of the controls in nCounter miRNA Expression assays.

Control Type	Number	Description	Use
Positive Controls	6	Probes that recognize synthetic mRNA targets included in the assay at specified concentrations (targets do not require ligation).	Positive controls used by the QC metrics in nSolver to confirm linear response to input amounts, and confirm that low input signal is above background.
Negative Controls	8	Probes that recognize synthetic mRNA targets not included in the assay (targets do not require ligation).	Negative controls used by the QC metrics in nSolver to determine background signal independent of ligation success. This can then be used to set threshold for defining expression of miRNA.
Ligation Positive Controls	3	Probes that recognize synthetic miRNA targets included in the Sample Preparation Kit (Sample Preparation Kit includes targets, bridges, and miRtags).	Ligation positive controls monitor ligation efficiency, independent of the miRNAs in the sample.
Ligation Negative Controls	3	Probes that recognize synthetic miRNA targets not included in the Sample Preparation Kit (no miRNA target).	Ligation negative controls monitor non-specific ligation.
mRNA Reference Controls	5	Probes that recognize endogenous mRNA targets commonly expressed in tissues.	Measurement of mRNA can be used to determine if there is cellular contamination in cell free samples. When miRNA is extracted from cells, even if ligation fails, the mRNA will be seen. It is not recommended that these mRNA be used as reference genes to normalize samples.
Spike-in Controls	5	Probes that recognize exogenous miRNA targets not included in the Sample Preparation Kit to monitor RNA isolation/purification steps upstream of the assay (optional, added by user if desired).	If user adds the appropriate small molecules to samples, they can be used to monitor steps of sample processing prior to the NanoString protocol.



and lymph node metastatic lesions from over 170 women. Their overall results helped define specific miRNA expression signatures that further characterize the complexity and diversity of TNBC as well as its metastasis<sup>35</sup>. Other clinical utility with the adoption of the NanoString nCounter platform lies in examining the hematopoeitic cell of origin of blood based miRNAs as well as the ability to explore expression profiling from individual blood cell types. Teruel-Montoya et al. used the nCounter Human miRNA expression assay to profile the miRNA content from RNA generated from platelets, T cells, B cells, granulocytes, and other hematological cell types. The authors suggested that such data could be used as a basis for interpretation of miRNA-disease association studies. For example, a miRNA elevated in acute myelogenous leukemia (AML), but absent or very low in normal granulocytes (as the authors demonstrated), might participate in the pathogenesis of AML. Further insight into the differential expression of miRNAs by blood cell type therefore is also relevant to understanding the systemic effects of blood cell delivered miRNAs<sup>36</sup>.

The comprehensive nCounter miRNA assays are also being used to discover biomarkers and develop expression signatures from biofluids including serum and plasma<sup>37-39</sup>. The discovery of these relatively stable, extracellular circulating miRNAs has generated a great deal of interest over conventional biomarkers, particularly their likelihood as candidates for non-invasive early disease monitoring and prognostic biomarkers of disease progression or resolution. To this end, Permuth-Wey J. et al. recently discovered novel miRNA biomarkers in plasma samples from patients with intraductal papillary mucinous neoplasms of the pancreas, a precursor to pancreatic ductal adenocarcinoma (PDAC), one of the most fatal cancers worldwide. These exciting developments using nCounter miRNA expression assays represent a potential new technology and methodology for a disease which, to date, lacks early detection or intervention at an operable stage<sup>40</sup>. Similarly, Armstrong et al. profiled miRNAs derived from matched FFPE tumors and plasma and urine-based biofluids from patients with bladder cancer using the nCounter miRNA expression assay. Bladder cancer remains one of the most expensive cancers to treat due to the high rate of local recurrence and the requirement for frequent long-term follow-up with invasive and uncomfortable cystoscopic evaluations.

There is a clear unmet need therefore for additional, improved urine- (or other bio-specimen) based alternatives to cystoscopy for the screening, initial evaluation, and follow-up of bladder cancer; miRNA profiling may afford such a novel methodology<sup>12</sup>. In other examples, Stylii *et al.* profiled miRNAs from cerebrospinal fluid (CSF) samples in patients with aneurysmal subarachnoid hemorrhage (SAH), prognosis of which is very poor relating to a high mortality rate. The findings indicated that temporal miRNA profiling using NanoString can detect differences between CSF from aneurysmal SAH and non-SAH patients. Moreover, the miRNA profile of CSF samples from patients who develop cerebral vasospasm; a common complication from patients suffering from SAH, may be distinguishable from those who do not<sup>41</sup>. Further, in a study of miRNA expression profiling in extracellular vesicles from serum samples, Rodosthenis et al. sought to uncover whether air pollution from particulate matter could be linked to cardiovascular disease and mortality through a deregulation of miRNAs. The group identified several miRNAs following exposure to particulate matter and interestingly conducted further studies and pathway analysis mapping to demonstrate the role of these identified miRNAs in cardiovascular disease related pathways<sup>42</sup>. The breadth of biological applications accommodated using the NanoString miRNA expression assays continues to expand. Recently, Wang et al. investigated the role of differentially-expressed miRNAs following human enterovirus 71 (EV71) infection, the causative agent of hand, foot, and mouth disease. The group adopted NanoString miRNA profiling to characterize the responses of serum miRNA profiles in normal samples as well as in samples from EV71 infected patients suffering symptoms ranging from mild to severe. 44 miRNAs were observed in patients with infections and were found to regulate a number of key genes involved in known signaling pathways associated with the diagnosis of infection<sup>43</sup>.

These are just a few examples of the nCounter miRNA profiling platform enabling identification and evaluation of differentiallyexpressed miRNA tissue and biofluid-based biomarkers.

In addition to the studies mentioned herein, NanoString maintains a comprehensive and frequently-updated online list of all of the publications that utilize NanoString assays, including miRNA, gene expression (mRNA), gene fusions, single cell and CNV. For more references and information, see NanoString's Publications page, <u>https://www.nanostring.com/scientific-content/publications.</u>



### **Summary and Conclusion**

Many diverse platforms have been employed to investigate miRNA profiles and signatures, from broad screening to targeted platforms and methods. NanoString's miRNA expression assays are a set of unique and powerful biomarker discovery and signature development tools, providing excellent specificity, low false-positive rates<sup>31</sup>, and a direct digital readout. With a simple workflow, these assays enable investigators to generate comprehensive expression data sets guickly with minimal hands-on manipulation. This white paper highlighted many of the diverse applications for NanoString miRNA assays, including cancer and hematopoiesis. nCounter miRNA profiling has enabled the identification of individual or multiple miRNA biomarkers that reflect changes in disease states or effects of therapeutic agents, demonstrating the versatility of the nCounter platform to both identify and validate differentially expressed miRNAs. Such studies are possible due to carefully curated assay designs that ensure each assay contains biologicallyrelevant content and appropriate controls.

# REFERENCES

- 1 Bartel DP (2004) MicroRNAs: Genomics, Biogenesis, Mechanism, and Function. *Cell* 116:281–297.
- Lee RC, Feinbaum RL, Ambros V (1993) The C. elegans heterochronic gene lin-4 encodes small RNAs with antisense complementarity to lin-14. *Cell* 75:843-854
- 3 Brennecke J, Hipfner DR, Stark A, Russell RB, Cohen SM (2003) bantam encodes a developmentally regulated microRNA that controls cell proliferation and regulates the proapoptotic gene hid in Drosophila. *Cell* 113:25–36.
- 4 Hwang H-W, Mendell JT (2006) MicroRNAs in cell proliferation, cell death, and tumorigenesis. Br J Cancer 94:776-780.
- 5 Calin GA, Croce CM (2006) MicroRNA signatures in human cancers. Nat Rev Cancer 6:857–866.
- 6 Kumar P, Dzso Z, MacKenzie C, Oestreicher J, Agoulnik S, Bryne M, Bernier F, Yanagimachi M, Aoshima K, Oda Y (2013) Circulating miRNA Biomarkers for Alzheimer's Disease. *PLoS* One 8(7):e69807.

- Wahid F, Khan T, Kim YY (2014) MicroRNA and diseases: Therapeutic potential as new generation of drugs. *Biochimie* 104:12–26.
- 8 Calin GA, Croce CM (2006) MicroRNA-cancer connection: The beginning of a new tale. *Cancer Res* 66:7390-7394.
- Deng S, Calin GA, Croce CM, Coukos G, Zhang L (2008) Mechanisms of microRNA deregulation in human cancer. *Cell Cycle* 7:2643–2646.
- **10** Schwarzenbach H, Nishida N, Calin G a, Pantel K (2014) Clinical relevance of circulating cell-free microRNAs in cancer. *Nat Rev Clin Oncol* 11:145–156.
- Beermann J, Piccoli MT, Viereck J, Thum T (2016) Non-coding RNAs in Development and Disease: Background, Mechanisms, and Therapeutic Approaches. *Physiological Reviews* 96 (4):1297-1325.
- 12 Armstrong DA, Green BB, Seigne JD, Schned AR, Marsit CJ (2015) MicroRNA molecular profiling from matched tumor and bio-fluids in bladder cancer. *Mol.Cancer* 14; 14:194.



- Geiss GK, Bumbarner RE, Birditt B, Dahl T, Dowidar N, Dunaway DL, Fell HP, Ferree S, Grogan T, James JJ, Maysuria M, Mitton JD, Oliveri P, Osborn JL, Peng T, Ratcliffe AL, Webster PJ, Davidson EH, Hood L, Dimitrov K (2008) Direct multiplexed measurement of gene expression with colorcoded probe pairs. *Nat Biotechnol* 26:317-325.
- Guerau-de-Arellano M, Alder H, Ozer HG., Lovett-Racke A, Racke MK (2012) miRNA profiling for biomarker discovery in Multiple Sclerosis: from microarray to deep sequencing. J. Neuroimmunol 248(1-2): 32-39.
- How C, Pintilie M, Bruce JP, Hui ABY, Clarke BA, Wong P, Yin S, Yan R, Waggott D, Boutros PC, Fyles A, Hedley DW, Hill RP, Milosevic M, Liu FF (2015) Developing a prognostic Micro-RNA signature for human cervical carcinoma. *PloS ONE* 10(4): e0123946.
- 16 Villaruz LC, Huang G, Romkes M, Kirkwood JM, Buch SC, Nukui T, Flaherty KT, Lee SJ, Wilson MA, Nathanson KL, Benos PV, Tawbi HA (2015) MicroRNA expression profiling predicts clinical outcome of carboplatin/paclitaxel-based therapy in metastatic melanoma treated on the ECOG-ACRIN trial E2603. *Clin Epigenetics* 7:58.
- Suh S, Yoo JY, Cui R, Kaur B, Huebner K, Lee TK, Ageilan RI, Croce CM (2014) FHIT suppresses epithelial-mesenchymal transition (EMT) and metastasis in lung cancer through modulation of microRNAs. *PloS Genet* 10(10):e1004652.
- **18** Murphy SJ, Lusardi TA, Phillips JI, Saugstad JA (2014) Sex differences in microRNA expression during developmentin rat cortex. *Neurochem Int* 77:24-32.
- 19 Karlsen TA, Jakobsen RB, Mikkelsen TS, Brinchmann JE (2014) microRNA-140 targets RALA and regulates chondrogenic differentiation of human mesenchymal stem cells by translational enhancement of SOX9 and ACAN. Stem *Cells Dev* 23:290–304.
- **20** Sardi SH, Glassner BJ, Chang JR, Yim SH (2014) Identification of critical biomarkers responsive to anti-autophagy therapies for pancreatic ductal adenocarcinoma through a performance analysis of miRNA platforms. *J Bioanal Biomed* S10:001.
- Ni B, Rajaram MVS, Lafuse WP, Landes MB, Schlesinger
  LS (2014) Mycobacterium tuberculosis decreases human
  macrophage IFN-γ responsiveness through miR-132 and miR-26a. J Immunol 193(9):4537-47.

- 22 Naqvi AR, Fordham JB, Khan A, Nares S (2013) MicroRNAs responsive to Aggregatibacter actinomycetemcomitans and Porphyromonas gingivalis LPS modulate expression of genes regulating innate immunity in human macrophages. *Innate Immun* 20(5):540-551.
- Yamamoto M, Singh A, Ruan J, Gauvreau GM, O'Byrne PM, Carlsten CR, FitzGerald JM, Boulet LP, Tebbutt SJ (2012)
   Decreased miR-192 expression in peripheral blood of asthmatic individuals undergoing an allergen inhalation challenge. *BMC Genomics* 13:655.
- 24 Fabbri M, Paone A, Calore F, Galli R, Gaudio E, Santhanam R, Lovat F, Fadda P, Mao C, Nuovo GJ, Zanesi N, Crawford M, Ozer GH, Wernicke D, Alder H, Caligiuri MA, Nana- Sinkam P, Perrotti D, Croce CM (2012) MicroRNAs bind to Toll-like receptors to induce prometastatic inflammatory response. *Proc Natl Acad Sci* 109(31):e2110-6.
- 25 Griffiths-Jones S, Grocock RJ, van Dongen S, Bateman A, Enright AJ (2006) miRBase: microRNA sequences, targets and gene nomenclature. *Nucleic Acids Res* 34:D140–D144.
- **26** Kozomara A, Griffiths-Jones S (2011) MiRBase: Integrating microRNA annotation and deep-sequencing data. *Nucleic Acids Res* 39:D152-157.
- Kozomara A, Griffiths-Jones S (2014) MiRBase: Annotating high confidence microRNAs using deep sequencing data. *Nucleic Acids Res* 42:D68-73.
- **28** Sundarbose K, Kartha R, Subramanian S (2013) MicroRNAs as biomarkers in cancer. *Diagnostics* 3:84–104.
- 29 Wahid F, Shehzad A, Khan T, Kim YY (2010) MicroRNAs: synthesis, mechanism, function, and recent clinical trials. *Biochim Biophys Acta* 1803:1231–1243.
- 30 Baker SC, Bauer SR, Beyer RP, Brenton JD, Bromley B, Burrill J, Causton H, Conley MP, Elespuru R, Fero M, Foy C, Fuscoe J, Gao X, Gerhold DL, Gilles P, Goodsaid F, Guo X, Hackett J, Hockett RD, Ikonomi P, Irizarry RA, Kawasaki ES, Kaysser-Kranich T, Kerr K, Kiser G, Koch WH, Lee KY, Liu C, Liu ZL, Lucas A, Manohar CF, Miyada G, Modrusan Z, Parkes H, Puri RK, Reid L, Ryder TB, Salit M, Samaha RR, Scherf U, Sendera TJ, Setterquist RA, Shi L, Shippy R, Soriano JV, Wagar EA, Warrington JA, Williams M, Wilmer F, Wilson M, Wolber PK, Wu X, Zadro R; External RNA Controls Consortium (2005) The External RNA Controls Consortium: a progress report. *Nat Methods* 2:731-734.



- 31 Chatterjee A, Leichter AL, Fan V, Tsai P, Purcell RV, Sullivan MJ, Eccles MR (2015) A cross comparison of technologies for the detection of microRNAs in clinical FFPE samples of hepatoblastoma patients. *Sci Rep.* 5:10438.
- 32 Valeri N, Braconi C, Gasparini P, Murgia C, Lampis A, Paulus-Hock V, Hart JR, Ueno L, Grivennikov SI, Lovat F, Paone A, Cascione L, Sumani KM, Veronese A, Fabbri M, Carasi S, Alder H, Lanza G, Gafa' R, Moyer MP, Ridgway RA, Cordero J, Nuovo GJ, Frankel WL, Rugge M, Fassan M, Groden J, Vogt PK, Karin M, Sansom OJ, Croce CM (2014) MicroRNA-135b promotes cancer progression by acting as a downstream effector of oncogenic pathways in colon cancer. *Cancer Cell* 25:469–483.
- 33 Hur K, Toiyama Y, Schetter AJ, Okugawa Y, Harris CC, Boland CR, Goel A (2015) Identification of a metastasis-specific microRNA signature in human colorectal cancer. *J Natl Cancer Inst* 107(3).
- 34 Saha P, Nanda R (2016) Concepts and targets in triplenegative breast cancer: recent results and clinical implications. *Ther Adv Med Oncol.* 8(5):351-9.
- 35 Cascione L, Gasparini P, Lovat F, Carasi S, Pulvirenti A, Ferro A, Alder H, He G, Vecchione A, Croce CM, Shapiro CL, Huebner K (2013) Integrated microRNA and mRNA signatures associated with survival in triple-negative breast cancer. PLoS One. 8(2):e55910.
- 36 Teruel-Montoya R, Kong X, Abraham S, Ma L, Kunapuli SP, Holinstat M, Shaw CA, McKenzie SE, Edelstein LC, Bray PF (2014) MicroRNA expression differences in human hematopoietic cell lineages enable regulated transgene expression. *PLoS One* 9(7):e102259.
- Jacob NK, Cooley JV, Yee TN, Jacob J, Alder H, Wickramasinghe P, Maclean KH, Chakravarti A (2013) Identification of sensitive serum microRNA biomarkers for radiation biodosimetry. *PLoS One* 8(2):e57603.

- Brown JN, Brewer HM, Nicora CD, Weitz KK, Morris MJ, Skabelund AJ, Adkins JN, Smith RD, Cho JH, Gelinas R (2014) Protein and microRNA biomarkers from lavage, urine, and serum in military personnel evaluated for dyspnea. *BMC Med Genomics* 7:58.
- 39 Pfeffer SR, Grossmann KF, Cassidy PB, Yang CH, Fan M, Kopelovich L, Leachman SA, Pfeffer LM (2015) Detection of Exosomal miRNAs in the Plasma of Melanoma Patients. *J Clin Med.* 4(12):2012-27.
- Permuth-Wey J, Chen DT, Fulp WJ, Yoder SJ, Zhang Y, Georgeades C, Husain K, Centeno BA, Magliocco AM, Coppola D, Malafa M (2015) Plasma MicroRNAs as Novel Biomarkers for Patients with Intraductal Papillary Mucinous Neoplasms of the Pancreas. Cancer Prev Res (Phila). 8(9):826-34.
- 41 Stylli SS, Adamides AA, Koldej RM, Luwor RB, Ritchie DS, Ziogas J, Kaye AH (2016) miRNA expression profiling of cerebrospinal fluid in patients with aneurysmal subarachnoid hemorrhage. *J Neurosurg.* Apr 29:1-9. [Epub ahead of print]
- **42** Rodosthenous RS, Coull BA, Lu Q, Vokonas PS, Schwartz JD, Baccarelli AA (2016) Ambient particulate matter and microRNAs in extracellular vesicles: a pilot study of older individuals. *Part Fibre Toxicol.* 8;13:13.
- **43** Wang RY, Weng KF, Huang YC, Chen CJ (2016) Elevated expression of circulating miR876-5p is a specific response to severe EV71 infections. *Sci Rep.* 6:24149.

# For more information, please visit nanostring.com

NanoString Technologies, Inc. 530 Fairview Avenue North Seattle, Washington 98109

T (888) 358-6266 nanostring.com F (206) 378-6288 info@nanostring.com Sales Contacts United States us.sales@nanostring.com EMEA: europe.sales@nanostring.com

Asia Pacific & Japan apac.sales@nanostring.com Other Regions info@nanostring.com

#### FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

© 2019 NanoString Technologies, Inc. All rights reserved. NanoString, NanoString Technologies, the NanoString logo and nCounter are trademarks or registered trademarks of NanoString Technologies, Inc., in the United States and/or other countries. All other trademarks and/or service marks not owned by NanoString that appear in this document are the property of their respective owners.

